English

## MERIT MAK<sup>TM</sup>NV INSTRUCTIONS FOR USE

### PRODUCT DESCRIPTION

The Merit MAK-NV (non-vascular) introducer system utilizes a small co-axial introducer with dilator and stiffening cannula and guide wire for placement of larger diameter guide wires in non-vascular procedures.

The MAK-NV consists of a combination of the following components, see unit label for specific product components. These components may be packaged in a single pouch or may be packaged separately.

- 6F Co-axial introducer with dilator with stiffening cannula
- 21 gauge introducer needle with trocar or chiba stylet
- Guide wire(s)

### INTENDED USE

The Merit MAK-NV is intended to facilitate introduction and placement of up to a 0.038” diameter guide wire for non-vascular procedures.

### WARNINGS

Withdrawal, pull back, or manipulation of the guide wire distal tip through the needle tip may result in breakage or embolization.

Do not advance the guide wire if resistance is met.

### CAUTIONS

- Device should be used for non-vascular procedures only.
- Consult accompanying documents.
- RX Only Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.
- This device is intended for single use only.
- Do not reuse or resterilize.
- Contents of unopened, undamaged package are sterile using ethylene oxide and non-pyrogenic.

### POTENTIAL COMPLICATIONS

The Merit MAK-NV should be used by persons knowledgeable of the risks involved and qualified in the procedures performed. Potential complications include risks normally associated with percutaneous diagnostic and/or interventional procedures.

Other complications include, but are not limited to:

- Edema
- Infection
- Hematoma
- Bleeding
- Hemothorax
- Hydrothorax
- Extravasation
- Thrombus formation
- Pseudo aneurysm formation
- Perforation or laceration of vessel or viscous
- Guide wire or catheter embolization
- Inflammation, necrosis or scarring
- Pain in region

### INSTRUCTIONS FOR USE

- Identify the insertion site and prepare the site using proper aseptic technique and local anesthesia as required.
- Insert the 21 gauge introducer needle using standard technique.
- Carefully advance the flexible end of the 0.018” (0.46mm) guide wire through needle. Advance the guide wire as far as appropriate. Verify correct positioning. WARNING: Do not advance guide wire if resistance is met.
- Remove needle while maintaining the 0.018” (0.46mm) guide wire in position. WARNING: Withdrawal, pull back, or manipulation of the guide wire distal tip through the needle tip may result in breakage or embolization. To avoid guide wire damage during manipulation, remove the introducer needle and proceed to step 5.
- Insert the co-axial introducer with dilator, with or without stiffening cannula, over the 0.018” (0.46mm) guide wire and advance to the desired position.
- Remove the dilator/stiffening cannula, leaving the introducer in position. The 0.018” (0.46mm) guide wire may be removed or left in position.
- Insert up to 0.038” (0.97mm) guide wire through introducer.
- Remove introducer leaving guide wire(s) in place.

### REUSE PRECAUTION STATEMENT

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/ or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/ or cause patient infection or crossinfection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

**CE** 0086  Manufacturer: Merit Medical Systems, Inc. South Jordan, Utah 84095 www.merit.com U.S.A. 1-801-253-1600 U.S.A. Customer Service 1-800-356-3748

**EC REP** Authorized Representative: Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland

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