TEMNO®

Coaxial Introducer Needle

EN	INSTRUCTIONS FOR USE	2
FR	MODE D'EMPLOI	3
IT	ISTRUZIONI PER L'USO	4
DE	GEBRAUCHSANWEISUNG	5
ES	INSTRUCCIONES DE USO	6
PT	INSTRUÇÕES DE UTILIZAÇÃO	7
PTBR	INSTRUÇÕES DE USO	8
NL	GEBRUIKSINSTRUCTIES	9
SV	BRUKSANVISNING	10
DA	BRUGSANVISNING	11
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	12
TR	KULLANMA TALİMATLARI	13
PL	INSTRUKCJA UŻYTKOWANIA	14
CS	NÁVOD K POUŽITÍ	15
BG	ИНСТРУКЦИИ ЗА УПОТРЕБА	16
HU	HASZNÁLATI ÚTMUTATÓ	17
RU	РУКОВОДСТВО ПО ЭКСПЛУАТАЦИИ	18
FI	KÄYTTÖOHJEET	19
ZH	使用说明	20
JA	使用説明	21


TEMNO®

Coaxial Introducer Needle

INSTRUCTIONS FOR USE

INDICATIONS FOR USE:

Temno® Coaxial Introducer Needles are intended to be used as needle guides in soft tissue biopsy procedures.

NOTE: Not intended for use in bone. Not compatible with MRI.

CAUTION: Use with imaging equipment such as ultrasound, X-ray, CT, etc.



Do not use if package is damaged.

DIRECTIONS FOR USE:

1. During the pre-procedure set-up, insert the biopsy device through the coaxial introducer needle to verify compatibility.
2. Inspect the needle for damage before use.
3. Adjust sliding skin stopper on coaxial introducer needle to desired calibration and lock in place.
4. Insert introducer needle until sliding stopper rests on skin surface. Image coaxial needle, if desired, to make sure the tip is at or near the edge of the region to be sampled.
5. Pull the stylet out of the cannula.
6. Insert a procedural biopsy needle into the introducer needle until the biopsy device touches the hub on the introducer needle. **NOTE:** Automated devices must be prepared to fire prior to insertion.
7. Perform the biopsy.
8. Withdraw the biopsy needle from the introducer needle and remove the specimen.
9. If another biopsy of the same tissue is required, repeat steps 6 and 7. **NOTE:** If changing the angle of the introducer needle is required, reinsert the stylet prior to repositioning the needle. Remove the stylet prior to taking the biopsy.
10. Remove the introducer needle from biopsy site.

CAUTION: For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or reesterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or reesterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Only clinicians who have been properly trained in this procedure, including patient preparation, sample preservation and the use of image guidance for percutaneous biopsy, should use this device.

SYMBOL	DESIGNATION
	Use By: YYYY-MM-DD
	Lot Number
	Catalog Number
	Sterilized Using Ethylene Oxide
	Do Not Use if Package is Damaged
	Single Use
Rx ONLY	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	For electronic copy scan QR code or go to www.merit.com/ifu and enter IFU ID Number. For printed copy, call U.S.A or E.U Customer Service.

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Merit Medical Systems, Inc.
South Jordan, Utah 84095 U.S.A.
U.S.A Customer Service 1-800-356-3748
Merit Medical Ireland Ltd,