INSTRUCTIONS FOR USE

Read Instructions Prior To Use

Rx ONLY: Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
Sterile (Ethylene Oxide). Non-pyrogenic.
For one time use only. Do not autoclave.

Indications

These Guiding Catheters are intended for use for intravascular introduction of interventional/diagnostic devices into coronary or peripheral vascular systems.

Contraindications

None known for guiding catheters.

Warnings

Discard guiding catheters after one procedure. Structural integrity and/or function may be impaired through reshaping, reuse, or resterilization. Use catheters prior to the expiration date on the package.

Precautions

- Store in a cool dark place.
- Do not use open or damaged packages.
- Do not resterilize the catheter.
- Do not autoclave. Exposure to temperatures above 54°C (130°F) may damage the catheter.
- Do not expose to solvents or ionizing radiation.
- Inspect the guiding catheter before use to verify that its size, shape, and condition are suitable for the specific procedure.
- Before use, flush all devices entering the blood vessel with sterile heparinized saline or a similar isotonic solution.
- If strong resistance is met during manipulation, discontinue the procedure and determine the cause of the resistance before proceeding. If the cause of the resistance cannot be determined, withdraw the catheter.
- Torquing the guiding catheter while kinked may cause damage which could result in possible separation along the catheter shaft. Should the guiding catheter shaft become kinked, withdraw the entire system (guiding catheter, guidewire, and catheter sheath introducer).
- Advancement, manipulation, and withdrawal of the guiding catheter should always be performed under fluoroscopic guidance.
- Extreme care must be taken to avoid damage to the vasculature through which the guiding catheter passes. The guiding catheter may occlude smaller vessels. Care must be taken to avoid complete blood flow blockage.
- Large internal lumen guiding catheters require less force on the syringe during injection.

Complications

Use of guiding catheters should be restricted to those specialists trained to perform procedures for which the product is indicated. Complications may occur at any time during or after the procedure and may include, but are not limited to, the following: hemorrhage, hematoma, allergic reaction, infection, embolism, blood vessel dissection or occlusion.

Recommended Procedure

1. Remove the guiding catheter from its packaging using sterile techniques.
2. Flush the guiding catheter lumen with a heparinized saline solution.
3. Appropriate anticoagulation and vasodilation therapy should be used.
4. Introduce the guiding catheter into the vasculature through the catheter sheath introducer, and/or over an indwelling guidewire using a percutaneous entry technique of choice.
5. Under fluoroscopic guidance, advance the guiding catheter over the guidewire or introducer until the desired position is attained.
6. Remove the guidewire prior to introduction of other intravascular devices or infusion of contrast agents.

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