


PRECISE NAVIGATION IN CHALLENGING ANATOMY

DELIVERY SYSTEMS


SwiftNINJA[®]

Steerable Microcatheter

The SwiftNINJA is the most advanced and only 180° articulating microcatheter on the market today. A straight tip catheter that articulates in opposing directions, the SwiftNINJA allows physicians to rapidly select vessels without the use of a guide wire.



Distal Tip Cross Section



Close-up of Steerable Tip



CATHETER COMPATIBILITY

French Size	Maximum Microsphere Size	Maximum Particle Size	Maximum Coil Size	Maximum Injection Pressure
2.9F - 2.4F	≤ 700 μm	≤ 710 μm	≤ 0.018" (0.46 mm)	1000 psi (6900 kPa)



"The frequent insertion and removal of guidewires during superselective catheterisation can be minimised by using a microcatheter with an articulating tip. Using a SwiftNINJA can potentially simplify and shorten the procedure!"

— Yasuaki Arai, MD, PhD

Chief, Department of Diagnostic Radiology, National Cancer Center of Japan

Catheter Preparation for Use

To activate the hydrophilic coating: Flush the catheter with heparinized saline via the standard luer BEFORE removal from the protective catheter sleeve.

To Unlock Steering Dial (white wheel): Hold the lower housing body and gently pull the steering dial toward the luer connector until a clicking sound is heard.



Prior to attempting catheter tip movement or deflection, you must engage the steering dial.

PRODUCT ORDERING INFORMATION

MERIT SWIFTNINJA STEERABLE MICROCATETER

Catalog Number	SwiftNINJA Proximal OD	SwiftNINJA Distal OD	SwiftNINJA Usable Length	SwiftNINJA Tip Shape	SwiftNINJA Hydrophilic Coating Length	
MIV-20500	2.9F	2.4F	125 cm	Articulating	80 cm	
	SwiftNINJA Distance From Tip to Distal Markerband	SwiftNINJA Distance Between 1st and 2nd Markerband	Minimum Guide Cath. ID	SwiftNINJA ID	Max. Guide Wire Size	SwiftNINJA Priming Vol.
	0.5 mm	13.5 mm	0.042" (1.07 mm)	0.021" (0.54 mm)	0.018" (0.46 mm)	0.49 mL

Distribuido por:



tecnología con sentido humano

ventas@cardioparts.com • www.cardioparts.com
Tels. (55) 9183-3933 (81) 8357-0779 (81) 8357-1903



1. Data on file

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733