




PERIPHERAL INTERVENTION

Merit PAK™

Pedal Access Kits



Your
Choice
for Limb
Salvage

Merit PAK™ Pedal Access Kits

Pedal access is a promising alternative access site for recanalization procedures in patients with Peripheral Artery Disease (PAD) and Critical Limb Ischemia (CLI)¹. Merit Medical has designed efficient kits to help physicians achieve Pedal Access: Introducing Merit PAK Pedal Access Kits.

Accurate Access

- **Echo-enhanced needles** for visibility under fluoroscopy and ultrasound
- **Clear hub needles** for visibility of flash

Mini Pedal Access Kit



Easy Insertion

- **Seamless transitions** from dilator to introducer for smooth insertions
- **Larger 4F (ID) PreludeEASE®** hydrophilic sheath for insertion of target treatment devices
- **Stiffened dilator** configuration available for tortuous anatomy

Merit PAK™ - a suite of pedal access kits including everything you need to access the arteries below the knees. The micro-needles, introducers, and dilators provide quick and easy access for limb salvage procedures.



The image displays the components of the Hydrophilic Pedal Access Kit against a white background with decorative water splashes in the top-left and bottom-right corners. The components include a long, thin micro-needle with a red and white handle, a red and white hemostasis valve adapter, a clear plastic introducer with a red stopcock, and a clear plastic dilator with a red stopcock. A clear plastic tube is also visible, connected to the micro-needle handle.

Hydrophilic Pedal Access Kit



Placement & Hemostasis

- **Super elastic Nitinol core** wire provides flexibility and kink resistance²
- **Hemostasis Valve Adapter** for management of injection, aspiration, and hemostasis

Merit PAK™ Pedal Access Kits

ORDERING INFORMATION

Catalog Number	Type	Introducer OD/ID	Sheath Working Length	Needle Gauge	Needle Length	Needle Hub & Type	Guide Wire Diameter	Guide Wire Length	Hemostasis Valve
PAKM4NPD4	Mini	4/3F	10cm	21	4cm	Clear/Echo	0.018"	40cm	Yes
PAKM4NPD7	Mini	4/3F	10cm	21	7cm	Clear/Echo	0.018"	40cm	Yes
PAKM5NPD4	Mini	5/3F	10cm	21	4cm	Clear/Echo	0.018"	40cm	Yes
PAKM5NPD7	Mini	5/3F	10cm	21	7cm	Clear/Echo	0.018"	40cm	Yes
PAKS4NPD4	Stiffened	4/3F	10cm	21	4cm	Clear/Echo	0.018"	40cm	Yes
PAKS4NPD7	Stiffened	4/3F	10cm	21	7cm	Clear/Echo	0.018"	40cm	Yes
PAKS5NPD4	Stiffened	5/3F	10cm	21	4cm	Clear/Echo	0.018"	40cm	Yes
PAKS5NPD7	Stiffened	5/3F	10cm	21	7cm	Clear/Echo	0.018"	40cm	Yes
PAKPHR47NPD4	Hydrophilic	6/4F	7cm	21	4cm	Clear/Echo	0.018"	40cm	Yes
PAKPHR47NPD7	Hydrophilic	6/4F	7cm	21	7cm	Clear/Echo	0.018"	40cm	Yes



MINI PEDAL ACCESS KIT

- Echo-enhanced, clear hub needle
- Introducer & dilator
- Nitinol/Palladium (NPD) Guide wire
- Hemostasis Valve Adapter

HYDROPHILIC PEDAL ACCESS KIT

- Echo-enhanced, clear hub needle
- Hydrophilic sheath introducer with hemostasis valve and a side arm
- Nitinol/Palladium (NPD) Guide wire

STIFFENED PEDAL ACCESS KIT

- Echo-enhanced, clear hub needle
- Introducer & stiffened dilator
- Nitinol/Palladium (NPD) Guide wire
- Hemostasis Valve Adapter

¹ El-Sayed H.F. Retrograde Pedal/Tibial Artery Access for Treatment of Infragenicular Arterial Occlusive Disease. Methodist DeBakey Cardiovasc J. 2013 Apr-Jun; 9(2): 73–78.

² Stockel D. Nitinol – Material with Unusual Properties. Endovascular Update, Issue 1. 1998

For a complete inventory please visit our website or talk to your Merit Medical Professional.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.208.4300
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Merit.com

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg
8002 25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Spain
+34 911238406

Sweden
020 792 445

Switzerland
(Deutsch)
+41 225180252
(Français)
+41 225948000
(Italiano)
+41 225180035
UK
0800 973 115