

DRAINAGE BAGS AND TUBING

Merit offers a selection of drainage bags and tubing to meet your specific drainage needs.

- Adjustable length tubing to fit a variety of patients
- Twist valve for easy disposal of fluids
- Internal anti-reflux valve prevents backflow
- Cloth backing for patient comfort
- Easy-to-read markings to identify fluid volume

Catalog Number	Description
DRAINAGE BAGS	
MDD600*	600 mL bag
MDD600PML*	600 mL bag with polypropylene luer
FZ624*	600 mL cloth-backed bag
K08-01907**	Jackson-Pratt 100 mL bulb
K10-03450**	MDD600 drainage bag with Jackson-Pratt® 100 mL bulb designed for vacuum drainage
ABSCESS DRAINAGE KIT	
K10-00230B†	OSHA-compliant closed waste system with 1,400 mL MDD bag, tubing and CRV Valve designed for quick and safe drainage elimination.
DRAINAGE CONNECTOR TUBING	
K10-04129‡	30 cm connecting tubing
K10-04130‡	42 cm connecting tubing
K10-04131‡	30 cm connecting tubing with one-way stopcock
K10-04132‡	42 cm connecting tubing with one-way stopcock
DRAINAGE TUBING	
FAS5248‡	48" large bore tubing with pinch clamp and non-vented spike
K08-02027‡	Tubing kit with barbed male luer
K10-04750†	48" large bore tubing with roller clamp and 16G x 1.5" (38.1 mm) needle
K08-02070†	96" large bore tubing with roller clamp and hose fitting for wall suction drainage

* Packaged 20 per box | ** Packaged 10 per box | † Packaged 15 per box | ‡ Packaged 25 per box
 Jackson-Pratt® is a registered trademark of Cardinal Health.



600 mL cloth-backed bag