

SPEED MEETS PRECISION

Embolization Gelatin Foam

Merit Medical embolization gelatin foam is uniformly cut and loaded into a syringe or torpedo cartridge for faster preparation and delivery.




S P E E D • P R E C I S I O N • C O N S I S T E N C Y™

PRODUCT INFORMATION

Ready-to-use products, just hydrate with 50:50 saline/contrast mixture:

EmboCube™ Embolization Gelatin



EmboCube™ Embolization Gelatin

- Uniformly cut gelatin foam cubes preloaded into Medallion® syringe

Torpedo™ Gelatin Foam



Torpedo™ Gelatin Foam

- Preshaped torpedo gelatin foam loaded into a cartridge with optional blunt stylet

EmboCube Embolization Gelatin Foam

CATALOG NUMBER	HYDRATED SIZE	WEIGHT	RECOMMENDED INJECTION SYRINGE	MINIMUM CATHETER ID	COMPATIBLE MERIT PRODUCT	COLOR CODE
EC2525	2.5 mm	25 mg	1 mL	0.024" (0.61 mm)	Merit Maestro® Microcatheter 2.8F - 2.8F Merit Maestro® Microcatheter 2.9F - 2.9F EmboCath® Plus Infusion Microcatheter	●
EC2550	2.5 mm	50 mg	1 mL	0.024" (0.61 mm)		●
EC25100	2.5 mm	100 mg	1 mL	0.024" (0.61 mm)		●
EC5025	5.0 mm	25 mg	3 mL	0.040" (1.02 mm)		●
EC5050	5.0 mm	50 mg	3 mL	0.040" (1.02 mm)		●
EC50100	5.0 mm	100 mg	3 mL	0.040" (1.02 mm)		●

Packaged 5/box.

Torpedo Gelatin Foam

CATALOG NUMBER	HYDRATED SIZE	LENGTH	RECOMMENDED INJECTION SYRINGE	MINIMUM CATHETER ID	COMPATIBLE MERIT PRODUCT	COLOR CODE
TOR2510	2.5 mm	10 mm	1 mL	0.027" (0.68 mm)	Merit Maestro® Microcatheter 2.9F - 2.9F EmboCath® Plus Infusion Microcatheter	●
TOR2520	2.5 mm	20 mm	1 mL	0.027" (0.68 mm)		●
TOR5010	5.0 mm	10 mm	3 mL	0.040" (1.02 mm)		●
TOR5020	5.0 mm	20 mm	3 mL	0.040" (1.02 mm)		●

Packaged 5/box.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
Main: +1.801.253.1600
Customer Service: +1.801.208.4300

Merit Medical Europe, Middle
East, & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

703936001_004 ID 072619