

Catalog #	Description	French	Length	Recom. Guide Wire	Side Ports	Bands	Band Spacing
<b>Performa® 2 Band</b>							
7742-2M65SH*	Pigtail Flush	4F	65 cm	.035" (0.89mm)	5	2	2 cm
7742-2M100SH*	Pigtail Flush	4F	100 cm	.035" (0.89 mm)	5	2	2 cm
7744-2M65SH*	Straight Flush	4F	65 cm	.035" (0.89 mm)	5	2	2 cm
7744-2M100SH*	Straight Flush	4F	100 cm	.035" (0.89 mm)	5	2	2 cm
7743-2M65SH*	Ultra Bolus Flush	4F	65 cm	.035" (0.89 mm)	5	2	2 cm
7743-2M100SH*	Ultra Bolus Flush	4F	100 cm	.035" (0.89 mm)	5	2	2 cm
7755-2M65	Mod. Hook Flush	4F	65 cm	.035" (0.89 mm)	6	2	2 cm
7755-2M100	Mod. Hook Flush	4F	100 cm	.035" (0.89 mm)	6	2	2 cm
7602-2M	Pigtail Flush	5F	100 cm	.035" (0.89 mm)	10	2	2 cm
7602-2M65	Pigtail Flush	5F	65 cm	.035" (0.89 mm)	10	2	2 cm
7674-2M	Ultra Bolus Flush	5F	100 cm	.035" (0.89 mm)	10	2	2 cm
7674-2M65	Ultra Bolus Flush	5F	65 cm	.035" (0.89 mm)	10	2	2 cm
7673-2M	Straight Flush	5F	100 cm	.035" (0.89 mm)	10	2	2 cm
7673-2M65	Straight Flush	5F	65 cm	.035" (0.89 mm)	10	2	2 cm
7603-2M100	Mod. Hook Flush	5F	100 cm	.035" (0.89 mm)	6	2	2 cm
7603-2M65	Mod. Hook Flush	5F	65 cm	.035" (0.89 mm)	6	2	2 cm
56535AVS2PIG**	Pigtail Flush	5F	65 cm	.035" (0.89 mm)	5	2	2 cm
510035AVS2PIG**	Pigtail Flush	5F	100 cm	.035" (0.89 mm)	5	2	2 cm
56535AVS2STR**	Straight Flush	5F	65 cm	.035" (0.89 mm)	5	2	2 cm
510035AVS2STR**	Straight Flush	5F	100 cm	.035" (0.89 mm)	5	2	2 cm
56535AVS2MHK**	Mod. Hook Flush	5F	65 cm	.035" (0.89 mm)	5	2	2 cm
510035AVS2MHK**	Mod. Hook Flush	5F	100 cm	.035" (0.89 mm)	5	2	2 cm

<b>Performa® 11 Band</b>							
7602-11M	Pigtail Flush	5F	100 cm	.035" (0.89 mm)	10	11	1 cm
7602-11M65	Pigtail Flush	5F	65 cm	.035" (0.89 mm)	10	11	1 cm
7673-11M100SH*	Straight Flush	5F	100 cm	.035" (0.89 mm)	10	11	1 cm
7673-11M65SH*	Straight Flush	5F	65 cm	.035" (0.89 mm)	10	11	1 cm
7603-11M100	Mod. Hook Flush	5F	100 cm	.035" (0.89 mm)	6	11	1 cm
7603-11M65	Mod. Hook Flush	5F	65 cm	.035" (0.89 mm)	6	11	1 cm

<b>Performa® 20 Band</b>							
7742-20M65SH*	Pigtail Flush	4F	65 cm	.035" (0.89 mm)	5	20	1 cm
7742-20M100SH*	Pigtail Flush	4F	100 cm	.035" (0.89 mm)	5	20	1 cm
7602-20M	Pigtail Flush	5F	100 cm	.035" (0.89 mm)	10	20	1 cm
7602-20M65	Pigtail Flush	5F	65 cm	.035" (0.89 mm)	10	20	1 cm
7602-20M100SH*	Pigtail Flush	5F	100 cm	.035" (0.89 mm)	5	20	1 cm
7602-20M65SH*	Pigtail Flush	5F	65 cm	.035" (0.89 mm)	5	20	1 cm
7673-20M100SH*	Straight Flush	5F	100 cm	.035" (0.89 mm)	10	20	1 cm
7673-20M65SH*	Straight Flush	5F	65 cm	.035" (0.89 mm)	10	20	1 cm
7603-20M100SH*	Mod. Hook Flush	5F	100 cm	.035" (0.89 mm)	6	20	1 cm
7603-20M65SH*	Mod. Hook Flush	5F	65 cm	.035" (0.89 mm)	6	20	1 cm
56535AVS20PIG**	Pigtail Flush	5F	65 cm	.035" (0.89 mm)	5	20	1 cm
510035AVS20PIG**	Pigtail Flush	5F	100 cm	.035" (0.89 mm)	5	20	1 cm
56535AVS20MHK**	Mod. Hook Flush	5F	65 cm	.035" (0.89 mm)	5	20	1 cm
510035AVS20MHK**	Mod. Hook Flush	5F	100 cm	.035" (0.89 mm)	5	20	1 cm

Catalog #	French	Length	Guide Wire	Recom. Guide Wire	Side Ports	Bands	Spacing
<b>Performa® Pediatric (11 mm pigtail curl)</b>							
7757-2M65	4F	65 cm	.035" (0.89 mm)	.041" (1.04 mm)	5	2	1 cm
7757-2M80	4F	80 cm	.035" (0.89 mm)	.041" (1.04 mm)	5	2	1 cm
7757-2M65SH*	4F	65 cm	.035" (0.89 mm)	.041" (1.04 mm)	5	2	1 cm
7757-2M80SH*	4F	80 cm	.035" (0.89 mm)	.041" (1.04 mm)	5	2	1 cm

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.253.1600  
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733

merit.com

Australia  
+61 (0) 3 8373 4943

Austria  
0800 295 374

Belgium  
0800 72 906 (Dutch)  
0800 73 172 (Français)

Brazil  
+55 11 4561 - 0062

Canada  
1.800.364.4370

China  
+86-10 8561 0788

Denmark  
80 88 00 24

Finland  
0800 770 586

France  
0800 91 60 30

Germany  
0800 182 0871

Hong Kong  
+852-2207 6888

India  
+91-80-41223376

Ireland (Republic)  
1800 553 163

Italy  
800 897 005

Luxembourg  
8002 25 22

Netherlands  
0800 022 81 84

Norway  
800 11629

Portugal  
308 801 034

Russia  
+7 495 221 89 02

Singapore  
+65 67536311

South Korea  
+8210 2205 0707

Spain  
+34 911238406

Sweden  
020 792 445

Switzerland  
800 11629  
+41 22 518 02 30  
(Français)  
+41 22 518 02 52

(Italiano)  
+41 22 518 00 35

UK  
0800 973 115

## RADIOPACITY

Catheter and marker bands are highly visible under fluoroscopy

## PLATINUM MARKER BANDS

These bands provide accurate sizing of vessel lumen before the interventional procedure. Marker band accuracy +/- .25 mm

## NYLON PEBAX SHAFT

Our kink resistant shaft provides optimal surface smoothness and thermal stability at body temperature.

## FLOW RATES<sup>1</sup>

4F		
Length	ml/sec	psi/kgf/cm <sup>2</sup>
65 cm	21	1192/84
100 cm	17	1190/84

5F		
Length	ml/sec	psi/kgf/cm <sup>2</sup>
65 cm	35	1192/84
90 cm	28	1142/80
100 cm	27	1196/84

4F Performa		
Length	ml/sec	psi/kgf/cm <sup>2</sup>
65 cm	21	1167/82
80 cm	15	1138/80

\* SH Configurations

\*\* AVS Configurations

<sup>1</sup> Flow rates achieved at 10.4cp contrast. Individual results may vary.