




# PUSH IT TO THE LIMIT™

*With Merit's Diagnostic Peripheral Catheters*



**Performa®**

**TORQUE  
PUSH  
FLOW**

# Performa<sup>®</sup>

---

Introducing Merit Medical's line of Diagnostic Peripheral Catheters. Designed to provide superior torque, increased shaft strength for better pushability and large inner lumens for improved flow rates.

Now you can ***PUSH IT TO THE LIMIT*** in today's endovascular environment with Merit's exceptional products.

## INNER-LUMEN SIZES

- 4F .042" (1.07mm)
- 5F .046" (1.17mm)
- 6F .054" (1.37mm)
- Yields higher flow rates at lower pressures
- Superior guide wire trackability

## STRAIN RELIEF

- Increased strength at catheter/hub interface
- Color coded for easy identification

## PRESSURE

- High pressure limits at given flow rates



## 4F CATHETER LINE

---

- ACCEPTS .038" (0.97MM) GUIDE WIRE FOR PHYSICIAN CONVENIENCE
- LARGE INNER DIAMETER .042" (1.07MM) PROVIDES EXCELLENT FLOW RATES
- BUMPER TIP ADDED FOR SOFTNESS AND PATIENT SAFETY

### WIRE BRAID DESIGN

- Greater torque for maximum control
- Improved longitudinal stability
- High degree of maneuverability
- Provides superior pushability

### WINGED POLYCARBONATE HUB

- Enhanced handling and control

### NYLON PEBAX SHAFT

- Kink Resistant
- Optimal surface smoothness
- Thermal stability at body temperature
- Thin wall design for improved flow rate

### RADIOPACITY


- Good shaft visualization under fluoroscopy
- Tip opacity ensures accurate placement and visualization

### TIP

- Soft tip for smooth atraumatic positioning
- Excellent tip memory
- Reforms accurately time after time

**HIGHER BRAID DENSITY  
FOR IMPROVED TORQUE  
SOFTER TIP  
INCREASED SHAFT RADIOPACITY  
PROVEN PERFORMER  
GREAT RADIOPACITY**

catalog number	French size	internal diameter	length	recommended guide wire	# of side ports	description
<b>HEADHUNTER</b>						
7715-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7715-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5575-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5575-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
7715-B0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7715-30	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5575-B3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5575-23	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
7715-C0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7715-50	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5575-C3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5575-33	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
7715-D0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7715-60	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5575-D3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5575-43	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
<b>NEWTON</b>						
7722-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7722-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5577-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5577-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
7722-B0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7722-20	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5577-B3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5577-23	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
7722-C0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7722-30	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5577-C3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5577-33	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
7722-D0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7722-40	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5577-D3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5577-43	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
7722-E0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7722-50	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5577-E3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5577-53	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
<b>BENTSON</b>						
7721-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7721-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5578-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5578-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
7721-B0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7721-20	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5578-B3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5578-23	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
7721-C0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7721-30	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5578-C3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5578-33	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	



catalog number	French size	internal diameter	length	recommended guide wire	# of side ports	description
----------------	-------------	-------------------	--------	------------------------	-----------------	-------------

**MANI**

7716-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7716-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5579-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5579-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	



MANI

**VERTEBRAL**


7717-A0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-	
7717-10	4	.042"/1.07mm	65cm	.038" (0.97mm)	-	bumper tip
7717-C0	4	.042"/1.07mm	80cm	.035" (0.89mm)	-	
7717-30	4	.042"/1.07mm	80cm	.038" (0.97mm)	-	bumper tip
7717-B0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7717-20	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5583-B3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-	
5583-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	-	
5583-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5583-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	



VERTEBRAL

**MODIFIED CEREBRAL**

7723-A0	4	.042"/1.07mm	90cm	.035" (0.89mm)	-	
7723-10	4	.042"/1.07mm	90cm	.038" (0.97mm)	-	bumper tip
5588-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5588-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	



MODIFIED CEREBRAL

**BERENSTEIN**

7718-C0	4	.042"/1.07mm	40cm	.035" (0.89mm)	-	
7718-30	4	.042"/1.07mm	40cm	.038" (0.97mm)	-	bumper tip
7718-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-	
7718-20	4	.042"/1.07mm	65cm	.038" (0.97mm)	-	bumper tip
7718-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7718-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5592-C3	5	.046"/1.17mm	40cm	.035" (0.89mm)	-	
5592-33	5	.046"/1.17mm	40cm	.038" (0.97mm)	-	
5592-B3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-	
5592-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	-	
5592-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5592-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	



BERENSTEIN

**STRAIGHT SELECTIVE**

7748-A0	4	.042"/1.07mm	110cm	.035" (0.89mm)	-	
7748-13	4	.042"/1.07mm	110cm	.038" (0.97mm)	-	bumper tip




STRAIGHT SELECTIVE

catalog number	French size	internal diameter	length	recommended guide wire	# of side ports	description
----------------	-------------	-------------------	--------	------------------------	-----------------	-------------

**MW2**

5590-13	5	.046"/1.17mm	40cm	.038" (0.97mm)	-	
5590-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	-	
5590-33	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	



**MODIFIED MW2**

7745-A0	4	.042"/1.07mm	40cm	.035" (0.89mm)	2	
7745-10	4	.042"/1.07mm	40cm	.038" (0.97mm)	2	bumper tip
7745-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	
7745-C0	4	.042"/1.07mm	100cm	.035" (0.89mm)	2	



**OSBORN**

7746-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7746-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
7746-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-	
7746-20	4	.042"/1.07mm	65cm	.038" (0.97mm)	-	bumper tip
7746-C0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7746-30	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip

1  
2



**HOOK 0.8**

7711-A0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	
7711-10	4	.042"/1.07mm	65cm	.038" (0.97mm)	2	bumper tip
5585-A3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2	
5585-13	5	.046"/1.17mm	65cm	.038" (0.97mm)	2	



**HOOK 1.0**

7711-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	
7711-20	4	.042"/1.07mm	65cm	.038" (0.97mm)	2	bumper tip
5585-B3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2	
5585-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	2	





**MODIFIED HOOK**

5585-C3	5	.046"/1.17mm	80cm	.035" (0.89mm)	-	
5585-D3	5	.046"/1.17mm	80cm	.035" (0.89mm)	-	
5585-E3	5	.046"/1.17mm	80cm	.035" (0.89mm)	-	

1  
2  
3



catalog number	French size	internal diameter	length	recommended guide wire	# of side ports	description	
<b>SIMMONS</b>							
7720-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	1	bumper tip
7720-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	-		
5576-F3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-		
5576-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-		
5576-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	-		
7720-B0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	2	bumper tip
7720-20	4	.042"/1.07mm	100cm	.038" (0.97mm)	-		
5576-G3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-		
5576-B3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-		
5576-23	5	.046"/1.17mm	100cm	.038" (0.97mm)	-		
7720-C0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	3	bumper tip
7720-30	4	.042"/1.07mm	100cm	.038" (0.97mm)	-		
5576-C3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-		
5576-33	5	.046"/1.17mm	100cm	.038" (0.97mm)	-		
7720-D0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	4	bumper tip
7720-40	4	.042"/1.07mm	100cm	.038" (0.97mm)	-		
5576-D3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-		
5576-43	5	.046"/1.17mm	100cm	.038" (0.97mm)	-		
<b>MODIFIED SIMMONS 1.5</b>							
7720-F0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-		
5576-E3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-		
<b>MODIFIED SIMMONS</b>							
7720-E0	4	.042"/1.07mm	80cm	.035" (0.89mm)	-		bumper tip
7720-50	4	.042"/1.07mm	80cm	.038" (0.97mm)	-		bumper tip
7720-70	4	.042"/1.07mm	110cm	.038" (0.97mm)	-		bumper tip
7720-80	4	.042"/1.07mm	120cm	.038" (0.97mm)	-		bumper tip
<b>COBRA 1</b>							
7710-A0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	1	bumper tip
7710-10	4	.042"/1.07mm	65cm	.038" (0.97mm)	2		
7710-D0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-		
7710-40	4	.042"/1.07mm	65cm	.038" (0.97mm)	-		
5584-A3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2		
5584-D3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-	2	bumper tip
5584-13	5	.046"/1.17mm	65cm	.038" (0.97mm)	-		
5584-43	5	.046"/1.17mm	65cm	.038" (0.97mm)	2		
7710-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2		
7710-E0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2		
7710-20	4	.042"/1.07mm	65cm	.038" (0.97mm)	-	2	bumper tip
7710-50	4	.042"/1.07mm	65cm	.038" (0.97mm)	-		
5584-B3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2		
5584-E3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-		
5584-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	-		
5584-53	5	.046"/1.17mm	65cm	.038" (0.97mm)	2	3	bumper tip
7710-C0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2		
7710-F0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-		
7710-30	4	.042"/1.07mm	65cm	.038" (0.97mm)	2		
7710-60	4	.042"/1.07mm	65cm	.038" (0.97mm)	-		
5584-C3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2	3	bumper tip
5584-F3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-		
5584-33	5	.046"/1.17mm	65cm	.038" (0.97mm)	-		
5584-63	5	.046"/1.17mm	65cm	.038" (0.97mm)	2		



catalog number	French size	internal diameter	length	recommended guide wire	# of side ports	description
----------------	-------------	-------------------	--------	------------------------	-----------------	-------------

**SHEPHERD HOOK 0.8**

7712-A0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	
7712-10	4	.042"/1.07mm	65cm	.038" (0.97mm)	2	bumper tip
5586-A3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2	
5586-13	5	.046"/1.17mm	65cm	.038" (0.97mm)	2	



SHEPHERD HOOK 0.8

**SHEPHERD HOOK 1.0**

7712-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	
7712-20	4	.042"/1.07mm	65cm	.038" (0.97mm)	2	bumper tip
5586-B3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2	
5586-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	2	



SHEPHERD HOOK 1.0

**SHEPHERD HOOK 1.0 MODIFIED**

7712-C0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	
7712-D0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-	
7712-30	4	.042"/1.07mm	65cm	.038" (0.97mm)	2	bumper tip
7712-40	4	.042"/1.07mm	65cm	.038" (0.97mm)	-	bumper tip
5586-43	5	.046"/1.17mm	65cm	.038" (0.97mm)	2	



SHEPHERD HOOK 1.0 MODIFIED

**RENAL DOUBLE CURVE**

7713-A0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	
7713-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-	
7713-10	4	.042"/1.07mm	65cm	.038" (0.97mm)	2	bumper tip
7713-20	4	.042"/1.07mm	65cm	.038" (0.97mm)	-	bumper tip
5587-A3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2	
5587-13	5	.046"/1.17mm	65cm	.038" (0.97mm)	-	
5587-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	2	



RENAL DOUBLE CURVE

**AMIR MOTARJEME**

7747-C0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-	
7747-D0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7747-E0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	
7747-F0	4	.042"/1.07mm	100cm	.035" (0.89mm)	2	
5591-A3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-	
5591-B3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5591-C3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2	
5591-D3	5	.046"/1.17mm	100cm	.035" (0.89mm)	2	



AMIR MOTARJEME

**AMIR MOTARJEME CANE**

7747-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-	
7747-20	4	.042"/1.07mm	65cm	.038" (0.97mm)	-	bumper tip
7747-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7747-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5591-E3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-	
5591-F3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	



AMIR MOTARJEME CANE



catalog number	French size	internal diameter	length	tip length	recommended guide wire	# of side ports	description
----------------	-------------	-------------------	--------	------------	------------------------	-----------------	-------------

### HOCKEY STICK

7752-10	4	.042"/1.07mm	100cm	4.5cm	.038" (0.97mm)	-	bumper tip
7543-13	5	.046"/1.17mm	100cm	4.5cm	.038" (0.97mm)	-	
7541-13	6	.054"/1.37mm	100cm	4.5cm	.038" (0.97mm)	-	




HOCKEY STICK

catalog number	French size	internal diameter	length	recommended guide wire	# of side ports	description
----------------	-------------	-------------------	--------	------------------------	-----------------	-------------

### REUTER

7714-A0	4	.042"/1.07mm	80cm	.035" (0.89mm)	2	
7714-10	4	.042"/1.07mm	80cm	.038" (0.97mm)	2	bumper tip
7714-B0	4	.042"/1.07mm	80cm	.035" (0.89mm)	-	
7714-20	4	.042"/1.07mm	80cm	.038" (0.97mm)	-	bumper tip
5589-A3	5	.046"/1.17mm	80cm	.035" (0.89mm)	-	
5589-B3	5	.046"/1.17mm	80cm	.035" (0.89mm)	2	
5589-13	5	.046"/1.17mm	80cm	.038" (0.97mm)	-	
5589-23	5	.046"/1.17mm	80cm	.038" (0.97mm)	2	



REUTER

### MIKAELSSON

7719-A0	4	.042"/1.07mm	80cm	.035" (0.89mm)	-	
7719-10	4	.042"/1.07mm	80cm	.038" (0.97mm)	-	bumper tip
5594-13	5	.046"/1.17mm	80cm	.038" (0.97mm)	-	



MIKAELSSON

### KA CATHETERS

5582-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	-	
5582-33	5	.052"/1.32mm	65cm	.038" (0.97mm)	-	high flow non-braided



KA CATHETERS

### KA2 CATHETERS

5582-13	5	.046"/1.17mm	40cm	.038" (0.97mm)	-	
5582-43	5	.046"/1.17mm	65cm	.038" (0.97mm)	-	
7514-13	6	.054"/1.37mm	40cm	.038" (0.97mm)	-	
7514-23	6	.054"/1.37mm	65cm	.038" (0.97mm)	-	



KA2 CATHETERS

### DVS A1 (10cm taper tip)

5595-13	5	.046"/1.12mm	65cm	.035" (0.89mm)	-	
5595-33	5	.046"/1.12mm	100cm	.035" (0.89mm)	-	
5595-23	5	.046"/1.12mm	65cm	.038" (0.97mm)	-	
5595-43	5	.046"/1.12mm	100cm	.038" (0.97mm)	-	



DVS A1

### DVS A2 (15cm taper tip)

5595-A3	5	.046"/1.12mm	65cm	.035" (0.89mm)	-	
5595-C3	5	.046"/1.12mm	100cm	.035" (0.89mm)	-	
5595-B3	5	.046"/1.12mm	65cm	.038" (0.97mm)	-	
5595-D3	5	.046"/1.12mm	100cm	.038" (0.97mm)	-	



DVS A2

catalog number	French size	internal diameter	length	recommended guide wire	# of side ports	description
<b>UHF SHEPHERD FLUSH</b>						
7602-A1	5	.052"/1.32mm	65cm	.035" (0.89mm)	6	non-braided
7602-B1	5	.052"/1.32mm	90cm	.035" (0.89mm)	6	non-braided
7602-C1	5	.052"/1.32mm	65cm	.038" (0.97mm)	6	non-braided
7602-D1	5	.052"/1.32mm	90cm	.038" (0.97mm)	6	non-braided



<b>ULTRA BOLUS FLUSH</b>						
7674-11	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	non-braided
7674-21	5	.052"/1.32mm	90cm	.035" (0.89mm)	10	non-braided
7674-31	5	.052"/1.32mm	65cm	.038" (0.97mm)	10	non-braided
7674-41	5	.052"/1.32mm	90cm	.038" (0.97mm)	10	non-braided



<b>PIGTAIL FLUSH</b>						
7749-A0	4	.042"/1.07mm	80cm	.035" (0.89mm)	10	braided
7749-A0	4	.042"/1.07mm	80cm	.038" (0.97mm)	10	braided



<b>ULTRA HIGH FLOW PIGTAIL</b>						
7602-41	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	non-braided
7602-51	5	.052"/1.32mm	90cm	.035" (0.89mm)	10	non-braided
7602-61	5	.052"/1.32mm	110cm	.035" (0.89mm)	10	non-braided
7602-71	5	.052"/1.32mm	65cm	.038" (0.97mm)	10	non-braided
7602-81	5	.052"/1.32mm	90cm	.038" (0.97mm)	10	non-braided
7602-91	5	.052"/1.32mm	110cm	.038" (0.97mm)	10	non-braided



<b>STRAIGHT FLUSH</b>						
7673-11	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	non-braided
7673-21	5	.052"/1.32mm	90cm	.035" (0.89mm)	10	non-braided
7673-31	5	.052"/1.32mm	65cm	.035" (0.89mm)	-	non-braided
7673-41	5	.052"/1.32mm	65cm	.038" (0.97mm)	10	non-braided
7673-51	5	.052"/1.32mm	90cm	.038" (0.97mm)	10	non-braided
7673-61	5	.052"/1.32mm	65cm	.038" (0.97mm)	-	non-braided




<b>MODIFIED HOOK FLUSH</b>						
7603-A1	5	.052"/1.32mm	65cm	.035" (0.89mm)	6	non-braided
7603-B1	5	.052"/1.32mm	90cm	.035" (0.89mm)	6	non-braided
7603-C1	5	.052"/1.32mm	65cm	.038" (0.97mm)	6	non-braided
7603-D1	5	.052"/1.32mm	90cm	.038" (0.97mm)	6	non-braided



catalog number	French size	internal diameter	length	recommended guide wire	# of side ports	bands/spacing	tip style
<b>PERFORMA VESSEL SIZING</b>							
<b>all side ports between tip and first band</b>							
7602-2M	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	2 bands 2cm spacing	Pigtail
7602-2M65	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	2 bands 2cm spacing	Pigtail
7674-2M	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	2 bands 2cm spacing	Ultra Bolus
7674-2M65	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	2 bands 2cm spacing	Ultra Bolus
7673-2M	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	2 bands 2cm spacing	Straight
7673-2M65	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	2 bands 2cm spacing	Straight
7603-2M100	5	.052"/1.32mm	100cm	.035" (0.89mm)	6	2 bands 2cm spacing	Mod. Hook
7603-2M65	5	.052"/1.32mm	65cm	.035" (0.89mm)	6	2 bands 2cm spacing	Mod. Hook
7742-2M65SH	4	.041"/1.04mm	65cm	.035" (0.89mm)	5	2 (2cm apart)	Pigtail
7742-2M100SH	4	.041"/1.04mm	100cm	.035" (0.89mm)	5	2 (2cm apart)	Pigtail
7744-2M65SH	4	.041"/1.04mm	65cm	.035" (0.89mm)	5	2 (2cm apart)	Straight
7744-2M100SH	4	.041"/1.04mm	100cm	.035" (0.89mm)	5	2 (2cm apart)	Straight
7743-2M65SH	4	.041"/1.04mm	65cm	.035" (0.89mm)	5	2 (2cm apart)	Ultra Bolus
7743-2M100SH	4	.041"/1.04mm	100cm	.035" (0.89mm)	5	2 (2cm apart)	Ultra Bolus
7755-2M65	4	.041"/1.04mm	65cm	.035" (0.89mm)	6	2 (2cm apart)	Mod. Hook
7755-2M100	4	.041"/1.04mm	100cm	.035" (0.89mm)	6	2 (2cm apart)	Mod. Hook

7602-11M	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	11 bands 1cm spacing	Pigtail
7602-11M65	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	11 bands 1cm spacing	Pigtail
7673-11M100SH	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	11 bands 1cm spacing*	Straight
7673-11M65SH	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	11 bands 1cm spacing*	Straight
7603-11M65	5	.052"/1.32mm	65cm	.035" (0.89mm)	6	11 bands 1cm spacing	Mod. Hook
7603-11M100	5	.052"/1.32mm	100cm	.035" (0.89mm)	6	11 bands 1cm spacing	Mod. Hook
7742-20M65SH	4	.041"/1.04mm	65cm	.035" (0.89mm)	5	20 bands 1cm spacing	Pigtail
7742-20M100SH	4	.041"/1.04mm	100cm	.035" (0.89mm)	5	20 bands 1cm spacing	Pigtail
7602-20M	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	20 bands 1cm spacing	Pigtail
7602-20M65	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	20 bands 1cm spacing	Pigtail
7602-20M65SH	5	.052"/1.32mm	65cm	.035" (0.89mm)	5	20 bands 1cm spacing*	Pigtail
7602-20M100SH	5	.052"/1.32mm	100cm	.035" (0.89mm)	5	20 bands 1cm spacing*	Pigtail
7673-20M65SH	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	20 bands 1cm spacing*	Straight
7673-20M100SH	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	20 bands 1cm spacing*	Straight
7603-20M100SH	5	.052"/1.32mm	100cm	.035" (0.89mm)	6	20 bands 1cm spacing*	Mod. Hook
7603-20M65SH	5	.052"/1.32mm	65cm	.035" (0.89mm)	6	20 bands 1cm spacing*	Mod. Hook



ULTRA BOLUS



PIGTAIL



STRAIGHT

MODIFIED HOOK

<b>PERFORMA VESSEL SIZING</b>							
<b>all side ports between bands</b>							
56535AVS2PIG	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	2 (2cm apart)	Pigtail
510035AVS2PIG	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	2 (2cm apart)	Pigtail
56535AVS2STR	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	2 (2cm apart)	Straight
510035AVS2STR	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	2 (2cm apart)	Straight
56535AVS2MHK	5	.052"/1.32mm	65cm	.035" (0.89mm)	6	2 (2cm apart)	Mod Hook
510035AVS2MHK	5	.052"/1.32mm	100cm	.035" (0.89mm)	6	2 (2cm apart)	Mod Hook
46535AVS20MHK	4	.041"/1.04mm	65cm	.035" (0.89mm)	6	20 (1cm apart)	Mod Hook
56535AVS20PIG	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	20 (1cm apart)	Pigtail
510035AVS20PIG	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	20 (1cm apart)	Pigtail
56535AVS20MHK	5	.052"/1.32mm	65cm	.035" (0.89mm)	6	20 (1cm apart)	Mod Hook
510035AVS20MHK	5	.052"/1.32mm	100cm	.035" (0.89mm)	6	20 (1cm apart)	Mod Hook



PIGTAIL



PIGTAIL

<b>SOFTOUCH PEDIATRIC VESSEL SIZING</b>							
4408-2M65	4	.041"/1.04mm	65cm	.035" (0.89mm)	5	2 bands (1cm apart)	Pigtail
4408-2M80	4	.041"/1.04mm	80cm	.035" (0.89mm)	5	2 bands (1cm apart)	Pigtail
4408-2M65SH	4	.041"/1.04mm	65cm	.035" (0.89mm)	5	2 bands (1cm apart)*	Pigtail
4408-2M80SH	4	.041"/1.04mm	80cm	.035" (0.89mm)	5	2 bands (1cm apart)*	Pigtail

\* Bands are 16mm closer to pigtail curve



Merit Medical Systems, Inc. • 1600 West Merit Parkway • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT

[www.merit.com](http://www.merit.com)

Merit Medical EUROPE, MIDDLE EAST AND AFRICA (EMEA) • Amerikalaan 42, 6199 AE Maastricht-Airport • The Netherlands • Tel: +31 43 358 82 40

Free phone for specific country:

Austria 0800 295374 • Belgium 0800 729 06 (Dutch) 0800 731 72 (French) • Denmark 80 88 00 24 • France 0800 91 60 30 • Germany 0800 182 08 71

Republic of Ireland 1800 553 163 • Luxembourg 8002 25 22 • Netherlands 0800 022 8184 • Sweden 020 79 24 45 • UK 0800 973 115 • Finland 00800 6374 8633

400459001/G ID 080811