

Performa[®]

DIAGNOSTIC CARDIOLOGY CATHETERS

MERITMEDICAL[®]

| Proven Performance

➤ PERFORMA® DIAGNOSTIC CARDIOLOGY CATHETERS

Introducing Merit Medical's Performa line of Diagnostic Cardiology Catheters, designed with all the torque you need, improved shaft strength for better pushability and large inner lumens for increased flow rates. The Performa line of catheters will help you reach your goal of excellent imaging in today's cardiology environment.

Performa®

NYLON MATERIAL

- Kink resistant
- Optimal surface smoothness
- Thermal stability at body temperature
- Thin wall design for improved flow rate

RADIOPACITY

- Good shaft visualization under fluoroscopy
- Tip opacity ensures accurate placement and visualization

WINGED POLYCARBONATE HUB

- Enhanced handling and control

WIRE BRAID DESIGN

- Greater torque for maximum control
- Improved longitudinal stability
- High degree of maneuverability
- Provides superior pushability

TIP DESIGNS

- Excellent tip memory
- Re-forms accurately time after time

BUMPER TIP

- Added softness to the 4F Performa catheter line for patient safety

FLAT WIRE BRAID

- Increased torque
- Greater stability
- Improved kink resistance



Proven Performer | Great Radiopacity | Improved Torque

4F PERFORMA CATHETER LINE

Accepts 0.038" (0.97 mm) guide wire for physician convenience

Bumper tip softness added for safe placement of tip into ostium

Large inner diameter 0.042" (1.07 mm) provides excellent flow rates

INNER LUMEN SIZES

- 4F 0.042" (1.07 mm)
- 5F 0.046" (1.17 mm)
- 5F 0.052" (1.32 mm) high flow
- 6F 0.054" (1.37 mm)
- 6F 0.059" (1.49 mm) high flow
- Yields higher flow rates at lower pressures
- Superior guide wire trackability

STRAIN RELIEF

- Increased strength at catheter/hub interface
- Color coded for easy identification

catalog number	French size	internal diameter	length	recommended guide wire	# side ports
----------------	-------------	-------------------	--------	------------------------	--------------

PEDIATRIC PIGTAIL CATHETERS - (FEP - FLUORINATED ETHYLENE PROPYLENE)

1310-21-0043	3	0.027"/.69 mm	40 cm	0.021" (0.53 mm)	4
1310-21-0053	3	0.027"/.69 mm	50 cm	0.021" (0.53 mm)	4
1310-21-0054	4	0.040"/1.02 mm	50 cm	0.021" (0.53 mm)	4
1310-21-0064	4	0.040"/1.02 mm	65 cm	0.021" (0.53 mm)	4
1310-25-0054	4	0.040"/1.02 mm	50 cm	0.025" (0.64 mm)	4
1310-25-0064	4	0.040"/1.02 mm	65 cm	0.025" (0.64 mm)	4
1310-25-0084	4	0.040"/1.02 mm	80 cm	0.025" (0.64 mm)	4
1310-25-0065	5	0.046"/1.17 mm	60 cm	0.025" (0.64 mm)	4
1310-25-0085	5	0.046"/1.17 mm	80 cm	0.025" (0.64 mm)	4
1308-25-0085	5	0.046"/1.17 mm	80 cm	0.025" (0.64 mm)	6
1308-35-0065	5	0.046"/1.17 mm	65 cm	0.035" (0.89 mm)	6
1308-35-0085	5	0.046"/1.17 mm	80 cm	0.035" (0.89 mm)	6
1308-35-0005	5	0.046"/1.17 mm	100 cm	0.035" (0.89 mm)	6
1308-35-0086	6	0.056"/1.42 mm	80 cm	0.035" (0.89 mm)	6
1308-35-0006	6	0.056"/1.42 mm	100 cm	0.035" (0.89 mm)	6
1308-38-0006	6	0.056"/1.42 mm	100 cm	0.038" (0.97 mm)	6
1308-35-0007	7	0.065"/1.65 mm	100 cm	0.035" (0.89 mm)	6
1308-38-0007	7	0.065"/1.65 mm	100 cm	0.038" (0.97 mm)	6





RADIAL APPROACH™

DIAGNOSTIC CARDIOLOGY CATHETERS

catalog number	French size	length	side ports
----------------	-------------	--------	------------

MIV™ VENTRICULOGAM CATHETER*

4RPS110	4	110 cm	6	braided pig
4RPS125	4	125 cm	6	braided pig
4RPL125	4	125 cm	6	braided pig
4RPL110	4	110 cm	6	braided pig
4RPS110-NB	4	110 cm	6	non-braided pig
4RPS125-NB	4	125 cm	6	non-braided pig
4RPL125-NB	4	125 cm	6	non-braided pig
4RPL110-NB	4	110 cm	6	non-braided pig
5RPS110	5	110 cm	6	braided pig
5RPS125	5	125 cm	6	braided pig
5RPL110	5	110 cm	6	braided pig
5RPL125	5	125 cm	6	braided pig
5RPS110-NB	5	110 cm	6	non-braided pig
5RPS125-NB	5	125 cm	6	non-braided pig
5RPL110-NB	5	110 cm	6	non-braided pig
5RPL125-NB	5	125 cm	6	non-braided pig
6RPS110	6	110 cm	6	braided pig
6RPS125	6	125 cm	6	braided pig
6RPL110	6	110 cm	6	braided pig
6RPL125	6	125 cm	6	braided pig
6RPS110-NB	6	110 cm	6	non-braided pig
6RPS125-NB	6	125 cm	6	non-braided pig
6RPL110-NB	6	110 cm	6	non-braided pig
6RPL125-NB	6	125 cm	6	non-braided pig



*Recommended guide wire is 0.038" (0.97 mm)

catalog number	French size	internal diameter	length	length of tip	# side ports	description
----------------	-------------	-------------------	--------	---------------	--------------	-------------

PERFORMA ULTIMATE™ RADIAL CATHETERS

ULTIMATE 1

410038ULT1	4	0.042" / 1.07 mm	100 cm	2.0 cm	-	bumper tip
4100381ULT1	4	0.042" / 1.07 mm	100 cm	2.0 cm	1	bumper tip
510038ULT1	5	0.046" / 1.17 mm	100 cm	2.0 cm	-	bumper tip
5100381ULT1	5	0.046" / 1.17 mm	100 cm	2.0 cm	1	bumper tip
5110381ULT1	5	0.046" / 1.17 mm	110 cm	2.0 cm	1	
610038ULT1	6	0.054" / 1.37 mm	100 cm	2.0 cm	-	
6100381ULT1	6	0.054" / 1.37 mm	100 cm	2.0 cm	1	
6110381ULT1	6	0.054" / 1.37 mm	110 cm	2.0 cm	1	



ULT1

ULTIMATE 2

410038ULT2	4	0.042" / 1.07 mm	100 cm	2.0 cm	-	bumper tip
4100381ULT2	4	0.042" / 1.07 mm	100 cm	2.0 cm	1	bumper tip
510038ULT2	5	0.046" / 1.17 mm	100 cm	2.0 cm	-	bumper tip
5100381ULT2	5	0.046" / 1.17 mm	100 cm	2.0 cm	1	bumper tip
5110381ULT2	5	0.046" / 1.17 mm	110 cm	2.0 cm	1	
610038ULT2	6	0.054" / 1.37 mm	100 cm	2.0 cm	-	
6100381ULT2	6	0.054" / 1.37 mm	100 cm	2.0 cm	1	
6110381ULT2	6	0.054" / 1.37 mm	110 cm	2.0 cm	1	



ULT2

ULTIMATE 3

410038ULT3	4	0.042" / 1.07 mm	100 cm	2.0 cm	-	bumper tip
4100382ULT3	4	0.042" / 1.07 mm	100 cm	2.0 cm	2	bumper tip
510038ULT3	5	0.046" / 1.17 mm	100 cm	2.0 cm	-	bumper tip
5100382ULT3	5	0.046" / 1.17 mm	100 cm	2.0 cm	2	bumper tip
610038ULT3	6	0.054" / 1.37 mm	100 cm	2.0 cm	-	
6100382ULT3	6	0.054" / 1.37 mm	100 cm	2.0 cm	2	
6110382ULT3	6	0.054" / 1.37 mm	110 cm	2.0 cm	2	



ULT3

ULTIMATE 4

(4 cm between primary & secondary curve)

410038ULT4-T40	4	0.042" / 1.07 mm	100 cm	2.0 cm	-	bumper tip
4100381ULT4-T40	4	0.042" / 1.07 mm	100 cm	2.0 cm	1	bumper tip
510038ULT4-T40	5	0.046" / 1.17 mm	100 cm	2.0 cm	-	bumper tip
5100381ULT4-T40	5	0.046" / 1.17 mm	100 cm	2.0 cm	1	bumper tip
5110381ULT4-T40	5	0.046" / 1.17 mm	110 cm	2.0 cm	1	
610038ULT4-T40	6	0.054" / 1.37 mm	100 cm	2.0 cm	-	
6100381ULT4-T40	6	0.054" / 1.37 mm	100 cm	2.0 cm	1	
6110381ULT4-T40	6	0.054" / 1.37 mm	110 cm	2.0 cm	1	



ULT4

(4.5 cm between primary & secondary curve)

410038ULT4-T45	4	0.042" / 1.07 mm	100 cm	2.0 cm	-	bumper tip
4100381ULT4-T45	4	0.042" / 1.07 mm	100 cm	2.0 cm	1	bumper tip
510038ULT4-T45	5	0.046" / 1.17 mm	100 cm	2.0 cm	0	bumper tip
5100381ULT4-T45	5	0.046" / 1.17 mm	100 cm	2.0 cm	1	bumper tip
5110381ULT4-T45	5	0.046" / 1.17 mm	110 cm	2.0 cm	1	
610038ULT4-T45	6	0.054" / 1.37 mm	100 cm	2.0 cm	0	
6100381ULT4-T45	6	0.054" / 1.37 mm	100 cm	2.0 cm	1	
6110381ULT4-T45	6	0.054" / 1.37 mm	110 cm	2.0 cm	1	

IMPRESS® MATERIAL

ULTIMATE 1

5IMP10038ULT1	5	0.046" / 1.17 mm	100 cm	2.0 cm	-	
5IMP100381ULT1	5	0.046" / 1.17 mm	100 cm	2.0 cm	1	



ULT1

Recommended guide wire is 0.038" (0.97 mm)



For ordering information call Merit Medical customer service 1-800-35-MERIT.



Merit Medical Systems, Inc. • 1600 West Merit Parkway • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT www.merit.com
Merit Medical EUROPE, MIDDLE EAST AND AFRICA (EMEA) • Amerikalaan 42, 6199 AE Maastricht-Airport • The Netherlands • Tel: +31 43 358 82 40
Free phone for specific country:
Austria 0800 295374 • Belgium 0800 729 06 (Dutch) 0800 731 72 (French) • Denmark 80 88 00 24 • France 0800 91 60 30 • Germany 0800 182 08 71
Republic of Ireland 1800 553 163 • Luxembourg 8002 25 22 • Netherlands 0800 022 8184 • Sweden 020 79 24 45 • UK 0800 973 115 • Finland 00800 6374 8633

ID 051711 400463001/F