

# Fluid Management

Options | Convenience | Confidence

**MERTMEDICAL**<sup>®</sup>

OPTIMIZE EFFICIENCY, REDUCE COST

# Fluid Management

## ➤ Integrate Fluid Management, Optimize Efficiency, and Reduce Cost

Merit Medical offers an integrated suite of fluid management products designed to effectively manage patient fluid levels, contrast media, and waste in diagnostic and interventional procedures. From tubing to transducers, fluid spikes to drainage bags, Merit Medical provides high quality angiographic accessories you can rely on to deliver safe and efficient fluid management.

### Manifolds & Stopcocks

#### Merit® Manifolds

- 1 to 5 ports
- Right or left handed orientation
- 200 psi or 500 psi
- 6" extension tubing & rotating Luer
- Block body or half body styles
- TRAM® and TRAM-P™ integral transducers and manifolds
- Devos™ manifold with integrated check relief valve

#### Marquis® Series Stopcock

- Large, easy-grip handle for effortless grip and turn
- Pressure ratings: 200, 500 & 1050 psi
- Luer options include rotator, swivel nut, fixed and slip
- 1-way and 3-way styles
- Off handles

### Contrast Management

#### Miser® Management System

- Exclusive ball and burette design to prevent air from entering fluid lines
- Dual one-way check valves' microbial barriers block cross contamination
- Automatic refill during each contrast withdrawal for quick and easy preparation

#### CT Contrast Transfer Sets

- Large and small bore configurations
- Multiple spike options
- Halkey-Roberts® swabbable valve



**OPTIONS** with a broad selection of fluid management products.

**CONVENIENCE** having all the items you need delivered together, streamlining clinical set up time.

**CONFIDENCE** with fluid and contrast management solutions that may help limit exposure for both patients and health care workers.



**Miser®**  
Contrast  
Management  
System

**Merit®**  
Flow Control Switch

## Tubing

### Merit® Tubing

- High Pressure contrast injection tubing rated at 1200, 900 and 500 psi (83, 62 and 34 Bar)
- Pressure Monitoring Tubing
- Fluid Administration Sets

## Waste Management

### Merit Disposal Depot™ Waste Bag

This closed system fully contains fluid waste for easy disposal, meeting OSHA guidelines. The bags come in many configurations for connecting directly to a manifold or for convenient placement on the back table.

- Large and small bore tubing available
- Fully customizable
- With and without reflux valves

## Transducers & Accessories

### Meritrans® Disposable Pressure Transducer

With reliable precision, customized configurations and a clear fluid path, the Meritrans can be customized to fit your clinical needs. Available in:

- Hand-held
- Pole-mount set-ups
- Cable Lengths: 24" or 48"
- Assembly with or without bonded stopcocks
- Mounting plate and pole clamp accessories

**Merit Disposal  
Depot™ Waste Bag**  
Closed Waste System

# Custom Kits

As a leading manufacturer of high quality products for cardiology, radiology, and other medical specialties, Merit Medical can offer custom fluid management to meet your unique clinical requirements, reducing in-house assembly and set up of components. Choose from:

- Syringes
  - Coronary Control Syringes (CCS™)
  - Medallion®
  - SmartTip™
  - Inject10n™
- PAL™ Pen and Labels
- Tubing Types & Lengths
- Manifold Set-up
- Pressure Infusor Bags

You will enjoy the simplicity of having all the items you need delivered together, streamlining clinical set up time.



Let Merit Medical customize a fluid management kit or complete procedure pack to meet your specific requirements. Customer Service (800) 35-MERIT or visit our website, [www.merit.com](http://www.merit.com).



**Merit Medical Systems, Inc.**  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1-801-253-1600  
1-800-35-MERIT  
[www.merit.com](http://www.merit.com)

**Merit Medical Europe, Middle East,  
and Africa (EMEA)**  
Amerikalaan 42, 6199 AE Maastricht-Airport  
The Netherlands  
Tel: +31 43 358 82 22

**Merit Medical Ireland Ltd.**  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733

#### Free phone for specific country

Austria 0800 295 374 | Belgium 0800 72 906 (Dutch)  
0800 73 172 (French) | Denmark 80 88 00 24  
France 0800 91 60 30 | Finland 0800 770 586 |  
Germany 0800 182 0871 | Ireland (Republic) 1800 553 163  
Italy 800 897 005 | Luxembourg 8002 25 22  
Netherlands 0800 022 81 84 | Norway 800 11629  
Sweden 020 792 445 | UK 0800 973 115