



HeartSpan[®]

Transseptal Needle



Product Features

Exceptionally responsive needle backed up with enhanced stiffness

Transparent handle allows direct visualization of bubbles before they travel distally

Ergonomic design offers excellent handling characteristics

Locking stylette provides guidance to the needle


HeartSpan
Transseptal
Needle Kit



(1) Transseptal Needle
(1) Stylette

HEARTSPAN Transseptal Needles

Part Number	Length	Curve Angle	Tip Gauge
FND-019-00	56 cm	50° Curve	22
FND-019-01	71 cm	50° Curve	21
FND-019-02	71 cm	86° Curve	21
FND-019-03	89 cm	50° Curve	21
FND-019-04	89 cm	86° Curve	21
FND-019-05	98 cm	50° Curve	21
FND-019-06	98 cm	86° Curve	21



Distribuido por:



CardioParts
Tecnología con sentido humano

ventas@cardioparts.com • www.cardioparts.com
Tels. (55) 9183-3933 (81) 8357-0779 (81) 8357-1903

