



PERIPHERAL INTERVENTION

MAK-NV™

**MAK-NV™**  
Introducer System

ACCESS

**Ideal for:**



Abscess Drainage



Biliary Drainage



Nephrostomy





Thoracentesis



Biliary Drainage




Nephrostomy



Paracentesis



Abscess Drainage



**Valve Adapter**  
for MAK-NV  
Introducer System

## START WITH QUALITY

Our MAK-NV (Mini Access Kit, Non-Vascular) Introducer System provides quick access and easy **VISUALISATION** into the drainage area with an echo-enhanced needle and radiopaque marker tip on the introducer. Superior tip design **DECREASES BUCKLING** and introducer peel-back for smooth insertion and **EXCELLENT TRANSITIONS**. The Valve Adapter is designed to be back-loaded over a guide wire and attached to the MAK-NV Introducer System to facilitate aspiration or infusion.



Introducer Needle  
with Chiba Stylet




Introducer Needle  
with Trocar Stylet

# MAK-NV™ (MINI ACCESS KIT, NON-VASCULAR) INTRODUCER SYSTEMS


The MAK-NV introducer system is designed for easy visualisation and quick access into the drainage area. For enhanced visibility, the device features an echo-enhanced needle and radiopaque marker tip on the introducer. Excellent transitions decrease buckling for a smooth insertion. Packaged 5 or 10 units per box.




Coaxial Introducer with Dilator and Stiffening Cannula



.018" (0.46 mm) Guide Wire with Platinum Tip



Introducer Needle with Trocar Stylet



Introducer Needle with Chiba Stylet



0.038" (0.97 mm) 3 mm J Guide Wire Heavy Duty Double Ended

Catalogue Number	Introducer	Dilator	Guide Wire 0.018" (0.46 mm) x 60 cm	Guide Wire Heavy Duty/Double Ended 0.038" (0.97 mm) 3mm J	Needle
MAK-NV-001	6F x 20 cm	4F	SS/PT	150 cm	21G x 15 cm Trocar
MAK-NV-002	6F x 20 cm	4F	NT/PT	150 cm	21G x 15 cm Trocar
MAK-NV-003	6F x 20 cm	4F	SS/PT	–	21G x 15 cm Trocar
MAK-NV-004	6F x 20 cm	4F	NT/PT	–	21G x 15 cm Trocar
MAK-NV-005*	–	–	–	–	21G x 15 cm Trocar
MAK-NV-006*	6F x 20 cm	4F	–	–	21G x 15 cm Trocar
MAK-NV-007*	6F x 20 cm	4F	–	–	–
MAK-NV-008	6F x 20 cm	4F	SS/PT	80 cm	21G x 15 cm Trocar
MAK-NV-009	6F x 20 cm	4F	NT/PT	80 cm	21G x 15 cm Trocar
MAK-NV-010	6F x 20 cm	4F	SS/PT	150 cm	21G x 15 cm Chiba
MAK-NV-011	6F x 20 cm	4F	NT/PT	150 cm	21G x 15 cm Chiba
MAK-NV-012	6F x 20 cm	4F	SS/PT	–	21G x 15 cm Chiba
MAK-NV-013	6F x 20 cm	4F	NT/PT	–	21G x 15 cm Chiba
MAK-NV-014*	–	–	–	–	21G x 15 cm Chiba
MAK-NV-015*	6F x 20 cm	4F	–	–	21G x 15 cm Chiba
MAK-NV-016	6F x 20 cm	4F	SS/PT	80 cm	21G x 15 cm Chiba
MAK-NV-017	6F x 20 cm	4F	NT/PT	80 cm	21G x 15 cm Chiba
MAK-NV-018	6F x 20 cm	4F	SS/PT	150 cm	21G x 20 cm Trocar
MAK-NV-019	6F x 20 cm	4F	NT/PT	150 cm	21G x 20 cm Trocar
MAK-NV-020*	–	–	–	–	21G x 20 cm Trocar
MAK-NV-021*	–	–	–	–	21G x 20 cm Chiba

SS = Stainless Steel PT = Platinum Tip NT = Nitinol



Packaged 10 per box.

\*Packaged 5 per box.

## VALVE ADAPTER FOR MAK-NV™ INTRODUCER SYSTEM

Designed to be back-loaded over a guide wire and attached to the MAK-NV Introducer System to facilitate aspiration or infusion. Will accept devices up to 4F. Packaged 10 per box.

**Catalogue Number** VA-40





## OPTIMISE YOUR DRAINAGE PROCEDURE

Trays, kits and packs are **EFFICIENT**, timesaving and cost-effective. Merit Medical offers numerous pre-arranged options, or we will work with you to create an assortment that **UNIQUELY** meets your needs and provides **EVERYTHING** you require to perform the drainage procedures you do most. Contact a Merit Medical Representative to learn more.



Understand. Innovate. Deliver.™



European Corporate Headquarters, Maastricht, The Netherlands

[www.Merit.com](http://www.Merit.com)

Distribuido por:



ventas@cardioparts.com • [www.cardioparts.com](http://www.cardioparts.com)  
Tels. (55) 9183-3933 (81) 8357-0779 (81) 8357-1903

**Merit Medical Europe, Middle East & Africa**

Amerikalaan 42, 6199 AE Maastricht-Airport, Netherlands  
+31 43 358 82 22

**Merit Medical Ireland Ltd.**

Parkmore Business Park West, Galway, Ireland  
+35 39 170 37 33

**Merit Medical France**

136 Ave. Joseph Kessel, 78960 Voisins Le Bretonneux, France  
+33 (0)1 30 57 24 24

**Merit Medical Systems, Inc.**

1600 West Merit Parkway, South Jordan, UT, 84095 USA  
+1 801 253 1600

**Asia, Pacific Rim, Latin America, Caribbean & Canada**

+1 801 208 4438

**Austria**

0800 295 374

**Belgium**

0800 72 906 (Dutch)  
0800 73 172 (Français)

**Denmark**

80 88 00 24

**Finland**

0800 770 586

**France**

0800 91 60 30

**Germany**

0800 182 0871

**Ireland (Republic)**

1800 553 163

**Italy**

800 897 005

**Luxembourg**

8002 25 22

**Netherlands**

0800 022 81 84

**Norway**

800 11629

**Portugal**

308 801 034

**Russia**

+7 495 221 89 02

**Sweden**

020 792 445

**Switzerland**

+41 225180252  
(Deutsch)

+41 225948000

(Français)

+41 225180035

(Italiano)

**UK**

0800 973 115