



# OPTIONS MATTER

Hemostasis Valves  
& Accessories



# Hemostasis Valves & Accessories

**NOW AVAILABLE IN THE U.S.**



## PhD™ Hemostasis Valve with Sidearm Tubing

Dual-seal hemostasis valve that simplifies procedures with an easy one-handed push-and-hold mechanism, allowing insertion of interventional equipment while controlling bleedback. Valve rated up to 400 psi.

### SPECIFICATIONS

Inner lumen: 0.096" (2.44mm)/7.3F  
Extension tube length: 23cm (9")



## PhD™ Hemostasis Valve

Dual-seal hemostasis valve that simplifies procedures with an easy push-and-hold mechanism, freeing your other hand to deliver interventional equipment. Valve rated up to 400 psi.

### SPECIFICATIONS

Inner lumen: 0.096" (2.44mm)/7.3F



## FLO40XR™

The FLO40XR Hemostasis Valve offers a shorter body and extension tubing as an additional feature. Valve rated up to 400 psi.

### SPECIFICATIONS

Inner lumen: 0.096" (2.44mm)/7.3F  
Extension tube length: 15cm (6")



## FLO50™

The FLO50 Hemostasis Valve offers a shorter body maintaining a large bore inner diameter rated up to 400 psi.

### SPECIFICATIONS

Inner lumen: 0.118" (3.0mm)/9F



## FLO30™

Merit Medical's FLO30 Hemostasis Valve is a short-body adapter providing maximum working length when required.

### SPECIFICATIONS

Inner lumen: 0.118" (3mm)/9F



## Honor®

The easy push-pull mechanism of the Honor Hemostasis Valve allows single-handed manipulation to open and close the valve while introducing devices into the vasculature.

### SPECIFICATIONS

Inner lumen: 0.094" (2.4mm)/7.2F



# Access the right hemostasis options



## AccessPLUS™

The AccessPLUS large bore hemostasis valve has optimal sealing surface for blood control.

### SPECIFICATIONS

Inner lumen: 0.120" (3.1mm)/9F



## DoublePlay™

Two large bore hemostasis ports allow easier manipulation and control of interventional devices during "kissing balloon" techniques.

### SPECIFICATIONS

Main inner lumen: 0.120" (3.1mm)/9F



## Passage®

The Passage represents the first of Merit's standard line of hemostasis valves.

### SPECIFICATIONS

Inner lumen: 0.097" (2.5mm)/7.4F



## MBA™ - Merit's Bleedback Alternative™

A premier hemostasis valve that features two seals. A conventional "touhy" silicone hemostasis valve provides security that the interventional device will stay in place. The secondary slit seal delivers a virtually bloodless procedure.

### SPECIFICATIONS

Inner lumen: 0.118" (3mm)/9F



## Prelude Choice™ Universal Hemostasis Valve Adapter

With two separate valve openings, the operator can track and control multiple devices while maintaining a secure seal that stops leaks. Compatible with standard luer devices, such as guide catheters and guiding sheaths with standard luer connections:

Hemostasis valves and accessories are also available with any Merit inflation device. Kits are available in standard catalog numbers as well as custom kit configurations.

Visit [merit.com/healthtrust/inflation-devices](http://merit.com/healthtrust/inflation-devices) for more information on the following devices:

basixTAU™ (pictured)


basixCOMPAK™

basixTOUCH™

basixTOUCH40™

Blue Diamond™

DiamondTOUCH™



# Ordering Information



## Merit Accessory Pak

Merit provides multiple hemostasis valve and accessory combinations tailored to your interventional needs.

Catalog Number	Hemostasis valve	Units per box
----------------	------------------	---------------

### STANDALONE HEMOSTASIS VALVES

MAP150	Access plus, large bore	25
MAP220	Double Play, large bore	25
MAP300	Honor, small bore	25
MAP350	Honor, small bore with bonded 20cm tubing and a stopcock	25
MAP400	MBA, large bore	25
MAP450	MBA, large bore with bonded 25.4cm high pressure tubing	25
MAP800	PhD, small bore	25
MAP850	PhD, small bore with bonded 23cm high pressure tubing	20
MAP900	Standalone FLO50 Valve	25
FLO30	FLO30, large bore	25
FLO40XR	Standalone FLO40XR Valve	20
K12-00040	Passage, small bore with bonded 15.2cm high pressure tubing	25

Body material: Polycarbonate

Catalog Number	Hemostasis valve	Units per box
----------------	------------------	---------------

### HEMOSTASIS VALVE ADAPTER

UHVA08	Prelude Choice Universal Hemostasis Valve Adapter	10
--------	---	----

Catalog Number	Hemostasis valve	Units per box	Accessories		
			Metal wire introducer	Plastic wire introducer	Torque device
<b>HEMOSTASIS VALVE MAP KITS</b>					
MAP101	Passage, small bore	25	X		
MAP102	Passage, small bore	25	X		X
MAP111	Access 9, large bore	25	X		
MAP112	Access 9, large bore	25	X		X
MAP113	Access 9, large bore	25		X	
MAP151	Access plus, large bore	25	X		
MAP152	Access plus, large bore	25	X		X
MAP153	Access plus, large bore	25		X	
MAP154	Access plus, large bore	25		X	X
MAP221	Double Play, large bore	25	X		
MAP222	Double Play, large bore	25	X		X
MAP301	Honor, small bore	25	X		
MAP302	Honor, small bore	25	X		X
MAP351	Honor, small bore with bonded 20cm tubing and a stopcock	25	X		
MAP352	Honor, small bore with bonded 20cm tubing and a stopcock	25	X		X
MAP353	Honor, small bore with bonded 50cm tubing and a stopcock	25	X		X
MAP403	MBA, large bore	25	X		X
MAP453	MBA, with bonded 25.4cm high pressure tubing	25	X		X
MAP801	PhD small bore	25	X		
MAP802	PhD small bore	25	X		X
MAP803	PhD small bore	25		X	
MAP804	PhD small bore	25		X	X
MAP851	PhD small bore with bonded 23cm high pressure tubing	20	X		
MAP852	PhD small bore with bonded 23cm high pressure tubing	20	X		X
MAP853	PhD small bore with bonded 23cm high pressure tubing	20		X	
MAP854	PhD small bore with bonded 23cm high pressure tubing	20		X	X
MAP952	FLO50 Valve, Metal Guide Wire Insertion Tool, & Torque Device	25	X		X
K05	Custom Kits				

### ACCESSORIES

MAP500	Torque Device 0.018", salmon	25			X
MAP501	Torque Device 0.018", green	25			X
MAP550	Metal Wire Insertion Tool	25	X		
MAP600	Plastic Wire Insertion Tool	25		X	

Body material: Polycarbonate

\* For additional compatibility information, please refer to Merit Medical's Prelude Choice™ sizing compatibility examples flier.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.253.1600  
1.800.35.MERIT

Distribuido por:

Tecnología con sentido humano

ventas@cardioparts.com • www.cardioparts.com  
Tels. (56) 9183-3933 (81) 8357-0779 (81) 8357-1903