




PERIPHERAL INTERVENTION

ANGIOGRAPHY

ANGIOGRAPHY



DIAGNOSTIC GUIDE WIRES

InQwire® Diagnostic Guide Wires	13
InQwire Diagnostic Guide Wires with Ringmaster™	15

HYDROPHILIC GUIDE WIRES

Merit Laureate® Hydrophilic Guide Wires	15
Tenor® Steerable Guide Wires	17

GUIDE WIRE ACCESSORIES

H2O Torq™ Torque Devices	18
SeaDragon™ Torque Device	18

DIAGNOSTIC PERIPHERAL CATHETERS

<i>Catheter Directory</i>	20
<i>Catheter Shape Illustrations</i>	21-22
Impress® Catheters	23-32
Impress Hydrophilic Catheters	33-36
Performa® Catheters	37-50

FLUID MANAGEMENT & TUBING

High Pressure Flow Control Switch	51
CT Transfer Sets	51

MICROCATHETERS

Merit Maestro® Microcatheters	52
-------------------------------	----

SYRINGES

Medallion® Syringes	53
Medallion Sword Handle Syringes	54
Medallion Frosted Field Syringes & Pen	54
Custom Medallion Syringes	55
VacLok® Vacuum Pressure Syringes	56

INQWIRE® DIAGNOSTIC GUIDE WIRES

The InQwire guide wire is precoated with PTFE over the entire wire surface for smooth, consistent coverage that resists flaking. Less friction means increased trackability, lubricity, and durability. Its precise J-tip memory does not require reloading the J straightener and the tight guide wire tolerances ensure catheter and sheath compatibility. The InQwire is also available with heparin coating to minimise clot formation.

Catalogue Number	Description	Size	Length cm	Core	Heparin Coated	# Per Box
STRAIGHT-TIP WIRES						
6676-A1	Straight	0.018" (0.46 mm)	150	F		10
6676-B3	Straight	0.018" (0.46 mm)	180	F		5
6676-C1	Straight	0.021" (0.53 mm)	150	F		10
6676-D3	Straight	0.021" (0.53 mm)	180	F		5
IQ35F150S	Straight	0.035" (0.89 mm)	150	F		10
IQH35F150S	Straight	0.035" (0.89 mm)	150	F	H**	10
IQ35F260S	Straight	0.035" (0.89 mm)	260	F		10
6690-41	Straight	0.035" (0.89 mm)	150	M		10
6690-61	Straight	0.038" (0.97 mm)	150	M		10
IQ38F150S	Straight	0.038" (0.97 mm)	150	F		10
IQ38F260S	Straight	0.038" (0.97 mm)	260	F		10
IQ35F150NLT	Newton; LT (13.5 cm)	0.035" (0.89 mm)	150	F		10
IQH35F150NLT	Newton; LT (13.5 cm)	0.035" (0.89 mm)	150	F	H**	10
IQ35F150NLLT	Newton; LLT (18.5 cm)	0.035" (0.89 mm)	150	F		10
IQH35F150NLLT	Newton; LLT (18.5 cm)	0.035" (0.89 mm)	150	F	H**	10
IQ35F150NLLLT	Newton; LLLT(23.5 cm)	0.035" (0.89 mm)	150	F		10
6677-B1	Newton; LT (13.5 cm) 3 mm J	0.035" (0.89 mm)	150	F		10
IQ35F80B	Bentson (23 cm)	0.035" (0.89 mm)	80	F		
IQ35F80BC	Bentson C	0.035" (0.89 mm)	80	F		10
IQ35F150B	Bentson (23 cm)	0.035" (0.89 mm)	150	F		10
IQH35F150B	Bentson (23 cm)	0.035" (0.89 mm)	150	F	H**	5/10*
IQ35F180B	Bentson (23 cm)	0.035" (0.89 mm)	180	F		10
IQH35F180B	Bentson (23 cm)	0.035" (0.89 mm)	180	F	H**	5
IQ35F260B	Bentson (23 cm)	0.035" (0.89 mm)	260	F		10
IQ38F150B	Bentson (23 cm)	0.038" (0.97 mm)	150	F		10
IQH38F150B	Bentson (23 cm)	0.038" (0.97 mm)	150	F	H**	10
IQ35F150BC	Bentson C (15 cm)	0.035" (0.89 mm)	150	F		10
IQ35F150BST	Bentson ST (10 cm)	0.035" (0.89 mm)	150	F		10
IQ35F180BC	Bentson C	0.035" (0.89 mm)	180	F		10
IQ35F180BST	Bentson ST (10 cm)	0.035" (0.89 mm)	180	F		10
IQ38F150BC	Bentson C	0.035" (0.89 mm)	150	F		10
IQ38F150BST	Bentson ST (10 cm)	0.038" (0.97 mm)	150	F		10
6680-61	Super Stiff	0.035" (0.89 mm)	80	F		10
6680-51	Super Stiff	0.035" (0.89 mm)	100	F		10
6680-11	Super Stiff	0.035" (0.89 mm)	150	F		10
6680-71	Super Stiff	0.035" (0.89 mm)	180	F		10
6680-23	Super Stiff	0.035" (0.89 mm)	260	F		5
6680-31	Super Stiff	0.038" (0.97 mm)	150	F		10
6680-43	Super Stiff	0.038" (0.97 mm)	180	F		5
ROSEN-TIP WIRES						
IQ35F80J1O5RS	1.5mm J (R)	0.035" (0.89 mm)	80	F		10
IQ35F150J1O5RS	1.5mm J (R)	0.035" (0.89 mm)	150	F		10
IQ35F180J1O5RS*	1.5mm J (R)	0.035" (0.89 mm)	180	F		5/10*
6609-33	1.5mm J (R)	0.035" (0.89 mm)	260	F		5
IQ38F150J1O5RS	1.5mm J (R)	0.038" (0.97 mm)	150	F		10

All "J" tip fixed core wires are finger-straightenable. FB = Firm Body, LT = Long Taper, ST = Short Taper, H = Heparin Coating, R = Rosen, F = Fixed Core, M = Moveable Core, TC = Tapered Core

*New guide wire configurations packaged 10 per box. Original guide wires packaged 5 per box.

**Heparin coated wires do not have CE Marking.

More options available through Merit's Cardiac Intervention catalogue. Ask for more details.

Catalogue Number	Description	Size	Length cm	Core	Heparin Coated	# Per Box
J-TIP WIRES						
IQ35F150J1O5	1.5 mm J	0.035" (0.89 mm)	150	F		10
IQH35F150J1O5	1.5 mm J	0.035" (0.89 mm)	150	F	H**	10
6693-A1	1.5 mm J	0.035" (0.89 mm)	150	M		10
IQ35F180J1O5	1.5 mm J	0.035" (0.89 mm)	180	F		10
IQ35F210J1O5	1.5 mm J	0.035" (0.89 mm)	210	F		10
IQ35F260J1O5	1.5 mm J Standard	0.035" (0.89 mm)	260	F		10
IQ35F260J1O5F	1.5 mm J FB	0.035" (0.89 mm)	260	F		10
6682-K1	3 mm J	0.018" (0.46 mm)	150	F		10
6682-L3	3 mm J	0.018" (0.46 mm)	180	F		5
6682-M1	3 mm J	0.021" (0.53 mm)	150	F		10
6682-N3	3 mm J	0.021" (0.53 mm)	180	F		5
6682-E1	3 mm J	0.025" (0.64 mm)	150	F		10
IQ35F80J3	3 mm J	0.035" (0.89 mm)	80	F		10
IQH35F150J3	3 mm J	0.035" (0.89 mm)	150	F	H**	10
IQ35F180J3*	3 mm J	0.035" (0.89 mm)	180	F		5/10*
IQ35F260J3	3 mm J	0.035" (0.89 mm)	260	F		10
6693-41	3 mm J	0.035" (0.89 mm)	150	M		10
6693-71	3 mm J	0.035" (0.89 mm)	150	M	H**	10
6647-11	3 mm J 7 cm TC	0.035" (0.89 mm)	145	M		10
IQ38F50J3	3 mm J	0.035" (0.89 mm)	50	F		10
IQ38F80J3	3 mm J	0.035" (0.89 mm)	80	F		10
IQ38F150J3	3 mm J	0.038" (0.97 mm)	150	F		10
IQH38F150J3	3 mm J	0.038" (0.97 mm)	150	F	H**	10
IQ38F260J3	3 mm J	0.038" (0.97 mm)	260	F		10
6693-61	3 mm J	0.038" (0.97 mm)	150	M		10
6684-61	3 mm J TC	0.038" (0.97 mm)	150	M		10
IQ35F150J3F	3 mm J FB	0.035" (0.89 mm)	150	F		10
IQH35F150J3F	3 mm J FB	0.035" (0.89 mm)	150	F	H**	10
IQ35F180J3F	3 mm J FB	0.035" (0.89 mm)	180	F		10
IQH35F180J3F	3 mm J FB	0.035" (0.89 mm)	180	F	H**	5
IQ38F150J3F	3 mm J FB	0.038" (0.97 mm)	150	F		10
IQ38F150J3S	3 mm J	0.038" (0.97 mm)	150	F		10
IQH38F150J3F	3 mm J FB	0.038" (0.97 mm)	150	F	H**	10
IQ38F180J3F	3 mm J FB	0.038" (0.97 mm)	180	F		10
IQ38F260J3F	3 mm J FB	0.038" (0.97 mm)	260	F		10
IQ35F80J6	6 mm J	0.035" (0.89 mm)	80	F		10
IQ35F150J6	6 mm J	0.035" (0.89 mm)	150	F		10
IQ35F180J6F	6 mm J FB	0.035" (0.89 mm)	180	F		10
IQH35F70J6	6 mm J	0.035" (0.89 mm)	70	F	H**	10
6694-21	6 mm J	0.035" (0.89 mm)	150	M		10
IQ38F150J6	6 mm J	0.038" (0.97 mm)	150	F		10
IQH38F70J6	6 mm J	0.038" (0.97 mm)	70	F	H**	10
IQ35100J15	15 mm J	0.035" (0.89 mm)	100	F		10
IQ35F150J15	15 mm J	0.035" (0.89 mm)	150	F		10
IQ35F180J15	15 mm J	0.035" (0.89 mm)	180	F		10
IQH35F150J15	15 mm J	0.035" (0.89 mm)	150	F	H**	10
IQ38F150J15	15 mm J	0.038" (0.97 mm)	150	F		10
8820-B1	15 mm J	0.035" (0.89 mm)	150	M		10
FIXED CORE, DOUBLE-ENDED						
IQ35F80J3SHD	3 mm J/Straight/Heavy Duty	0.035" (0.89 mm)	80	F		10
IQ35F150J105S	1.5 mm J/Straight	0.035" (0.89 mm)	150	F		10
IQ35F150J3S	3 mm J/Straight	0.035" (0.89 mm)	150	F		10
IQH35F150J3S	3 mm J/Straight	0.035" (0.89 mm)	150	F	H**	10
IQ35F150J3SHD	3 mm J/Straight/Heavy Duty	0.035" (0.89 mm)	150	F		10
IQ35F180J3S	3 mm J/Straight/Double End	0.035" (0.89 mm)	180	F		10
IQ38F150J3S	3 mm J	0.038" (0.97 mm)	150	F		10

All "J" tip fixed core wires are finger-straightenable. FB = Firm Body, LT = Long Taper, ST = Short Taper, H = Heparin Coating, R = Rosen, F = Fixed Core, M = Moveable Core, TC = Tapered Core

*New guide wire configurations packaged 10 per box. Original guide wires packaged 5 per box.

**Heparin coated wires do not have CE Marking.

More options available through Merit's Cardiac Intervention catalogue. Ask for more details.

INQWIRE® DIAGNOSTIC GUIDE WIRES WITH RINGMASTER™

Catalogue Number	Description	Size	Length cm	Core	# Per Box
------------------	-------------	------	-----------	------	-----------

STRAIGHT-TIP WIRES

IQR35F150B	Bentson (23 cm)	0.035" (0.89 mm)	150	F	10
IQR35F150BST	Bentson Short Taper (10 cm)	0.035" (0.89 mm)	150	F	10

J - TIP WIRES

IQR35F150J3	J 3 mm	0.035" (0.89 mm)	150	F	10
IQR35F260J3	J 3 mm	0.035" (0.89 mm)	260	F	10
IQR35F150J3F	J 3 mm/Firm	0.035" (0.89 mm)	150	F	10
IQR38F150J3	J 3 mm	0.035" (0.97 mm)	150	F	10
IQR38F150J3F	J 3 mm/Firm	0.035" (0.97 mm)	150	F	10

FIXED CORE: DOUBLE ENDED

IQR35F150J3S	J 3 mm/Straight	0.035" (0.89 mm)	150	F	10
--------------	-----------------	------------------	-----	---	----



Ringmaster™
Guide Wire Basin

Heparin coating not available with InQwire guide wires packaged in Ringmaster.

MERIT LAUREATE® HYDROPHILIC GUIDE WIRES

The Merit Laureate Hydrophilic Guide Wire facilitates the placement of devices during diagnostic and interventional procedures, navigating tortuous vessels. The hydrophilic coating offers enduring lubricity and consistent surface coating adhesion. The Merit Laureate was designed to reduce wire whipping, minimise friction, and to promote rapid catheter exchanges.

Catalogue Number	Description	Size	Length cm	Core	# Per Box
------------------	-------------	------	-----------	------	-----------

STANDARD WIRES

LWSTDA35150	Standard/Angled	0.035" (0.89 mm)	150	F	5
LWSTDA35180	Standard/Angled	0.035" (0.89 mm)	180	F	5
LWSTDA3880	Standard/Angled	0.038" (0.97 mm)	80	F	5
LWSTDA38120	Standard/Angled	0.038" (0.97 mm)	120	F	5
LWSTDA38150	Standard/Angled	0.038" (0.97 mm)	150	F	5
LWSTDA38180	Standard/Angled	0.038" (0.97 mm)	180	F	5
LWSTDS3580	Standard/Straight	0.035" (0.89 mm)	80	F	5
LWSTDS35150	Standard/Straight	0.035" (0.89 mm)	150	F	5
LWSTDS35180	Standard/Straight	0.035" (0.89 mm)	180	F	5
LWSTDS35220	Standard/Straight	0.035" (0.89 mm)	220	F	5
LWSTDS35260EX	Standard/Straight/Exchange	0.035" (0.89 mm)	260	F	5
LWSTDS38260EX	Standard/Straight/Exchange	0.038" (0.97 mm)	260	F	5
LWSTDS38150	Standard/Straight	0.038" (0.97 mm)	150	F	5
LWSTDS38180	Standard/Straight	0.038" (0.97 mm)	180	F	5
LWSTDA35260EX	Standard/Angled/Exchange	0.035" (0.89 mm)	260	F	5
LWSTDA38260EX	Standard/Angled/Exchange	0.038" (0.97 mm)	260	F	5

More options available through Merit's Cardiac Intervention catalogue. Ask for more details.

MERIT LAUREATE® HYDROPHILIC GUIDE WIRES

Catalogue Number	Description	Size	Length cm	Core	# Per Box
STIFF WIRES					
LWSTFA3580	Stiff/Angled	0.035" (0.89 mm)	80	F	5
LWSTFA35150	Stiff/Angled	0.035" (0.89 mm)	150	F	5
LWSTFA35180	Stiff/Angled	0.035" (0.89 mm)	180	F	5
LWSTFA35220	Stiff/Angled	0.035" (0.89 mm)	220	F	5
LWSTFA35260EX	Stiff/Angled/Exchange	0.035" (0.89 mm)	260	F	5
LWSTFA3880	Stiff/Angled	0.038" (0.97 mm)	80	F	5
LWSTFA38150	Stiff/Angled	0.038" (0.97 mm)	150	F	5
LWSTFA38180	Stiff/Angled	0.038" (0.97 mm)	180	F	5
LWSTFS3580	Stiff/Straight	0.035" (0.89 mm)	80	F	5
LWSTFS35150	Stiff/Straight	0.035" (0.89 mm)	150	F	5
LWSTFS35180	Stiff/Straight	0.035" (0.89 mm)	180	F	5
LWSTFS35220	Stiff/Straight	0.035" (0.89 mm)	220	F	5
LWSTFS35260EX	Stiff/Straight/Exchange	0.035" (0.89 mm)	260	F	5
LWSTFS38150	Stiff/Straight	0.038" (0.97 mm)	150	F	5
LWSTFS38180	Stiff/Straight	0.038" (0.97 mm)	180	F	5



The Merit Laureate® Hydrophilic Guide Wire has an improved hydrophilic coating that offers enduring lubricity and consistent turn-for-turn responsiveness.

Catalogue Number	Description	Size	Length cm	Core	# Per Box
SMALL DIAMETER SPECIALTY WIRE					
LWSTDA1880	Standard/Angled	0.018" (0.46 mm)	80	F	5
LWSTDA18150	Standard/Angled	0.018" (0.46 mm)	150	F	5
LWSTDA18180	Standard/Angled	0.018" (0.46 mm)	180	F	5
LWSTDA18260EX	Standard/Angled/Exchange	0.018" (0.46 mm)	260	F	5
LWSTDS18150	Standard/Straight	0.018" (0.46 mm)	150	F	5
LWSTDS18180	Standard/Straight	0.018" (0.46 mm)	180	F	5
LWSTDS18260EX	Standard/Straight/Exchange	0.018" (0.46 mm)	260	F	5
LWSTDA25150	Standard/Angled	0.025" (0.64 mm)	150	F	5
LWSTDA25180	Standard/Angled	0.025" (0.64 mm)	180	F	5
LWSTDA25260EX	Standard/Angled/Exchange	0.025" (0.64 mm)	260	F	5
LWSTDS25150	Standard/Straight	0.025" (0.64 mm)	150	F	5
LWSTDS25180	Standard/Straight	0.025" (0.64 mm)	180	F	5
LWSTDS25260EX	Standard/Straight/Exchange	0.025" (0.64 mm)	260	F	5


More options available through Merit's Cardiac Intervention catalogue. Ask for more details.

TENOR® STEERABLE GUIDE WIRES

The Tenor Steerable Guide Wire offers exceptional pushability and torqueability to facilitate precise navigation in challenging anatomy during embolic procedures. Packaged with an insertion tool and torque device.

Catalogue Number	Description	Size	Length cm	Tip Material	Core	# Per Box
TNR2411	Angled Tip	0.014" (0.36 mm)	165	Platinum Coiled	F	1
TNR2412	Straight Tip	0.014" (0.36 mm)	165	Platinum Coiled	F	1
TNR2811	Angled Tip	0.018" (0.46 mm)	165	Platinum Coiled	F	1
TNR2812	Straight Tip	0.018" (0.46 mm)	165	Platinum Coiled	F	1

Tenor Advanced Guide Wire Construction



TNR2411

H2O TORQ™ TORQUE DEVICE

The H2O Torq is an ergonomic torque device operable with one hand, designed to enhance manipulation of the guide wire tip. The single component plastic design can be used on hydrophilic or PTFE-coated wires and allows easy reposition and removal for rapid device exchanges. Packaged individually, 10 units per box.



TD100



TD200

Catalogue Number	Description
TD100	Orange device; Guide wires from 0.025" (0.64 mm) to 0.040" (1.02 mm)
TD200	Yellow device; Guide wires from 0.010" (0.25 mm) to 0.020" (0.51 mm)


SEADRAGON™ TORQUE DEVICE

The SeaDragon is an ergonomic torque device operable with one hand, designed to enhance manipulation of the guide wire tip. The single component plastic design can be used on hydrophilic or PTFE-coated wires and allows easy reposition and removal for rapid device exchanges. Packaged individually, 10 units per box.



Catalogue Number	Description
TD500	Guide wires from 0.018" (0.46 mm) to 0.038" (0.97 mm)



DIAGNOSTIC PERIPHERAL CATHETERS





DIAGNOSTIC PERIPHERAL CATHETERS

<i>Catheter Shape Illustrations</i>	21-22	KA2	34
		Mani	34
		Multipurpose	34
IMPRESS®		RBI	35
Headhunter	23	Simmons	35
Newton	23	Modified Hook	35
Bentson	24	Straight Selective	36
Mani	24	Vertebral	36
Vertebral	24		
Modified Cerebral	24	PERFORMA®	
Berenstein	25	Headhunter	37
MW2	25	Newton	38
Simmons	25	Bentson	39
Cobra	26	Mani	39
Hook	26	Vertebral	39
Shepherd Hook	27	Modified Cerebral	40
Renal Double Curve	27	Berenstein	40
Amir Motarjeme®	27	Straight Selective	40
Reuter	28	MW2	41
Mikaelsson	28	Osborn	41
KA2	28	Hook	41
Uterine Artery Catheter	28	Modified Hook Selective	42
Hockey Stick	28	Simmons	42
Modified Hook Selective	29	Cobra	43
RBI (Renal, Bifurcation, Iliac)	29	Shepherd Hook	44
RIM – Mini RIM Short-Tip	29	Renal Double Curve	44
Straight Selective	29	Hockey Stick	44
Osborn	30	Amir Motarjeme®	45
Multipurpose A1	30	Reuter	45
30 cm Catheters	30	Mikaelsson	46
Shepherd Flush	30	KA, KA2	46
Ultra Bolus Flush	31	DVS A1, DVS A2	46
Modified Bolus Flush	31	Shepherd Flush	47
Pigtail Flush	31	Ultra Bolus Flush	47
Straight Flush	32	Ultra High Flow Pigtail	47
Modified Hook Flush	32	Pigtail Flush	47
		Straight Flush	48
IMPRESS® HYDROPHILIC		Modified Hook Flush	48
Bentson	33	Vessel-Sizing, Adult	49-50
Berenstein	33		
Cobra	33		
Headhunter	34		


- a. HEADHUNTER 1
- b. HEADHUNTER 3
- c. HEADHUNTER 5
- d. HEADHUNTER 6





- a. SIMMONS 1
- b. SIMMONS 2
- c. SIMMONS 3
- d. SIMMONS 4
- e. MODIFIED SIMMONS 1.5
- f. MODIFIED SIMMONS





- a. NEWTON 1
- b. NEWTON 2
- c. NEWTON 3
- d. NEWTON 4
- e. NEWTON 5
- f. OSBORN 1
- g. OSBORN 2




- a. BENTSON 1
- b. BENTSON 2
- c. BENTSON 3
- d. MANI




- a. BERENSTEIN
- b. VERTEBRAL
- c. MODIFIED CEREBRAL
- d. MW2
- e. MODIFIED MW2
- f. DVS A1
- g. DVS A2
- h. HOCKEY STICK




- a. KA
- b. KA2
- c. HOOK 0.8
- d. HOOK 1.0
- e. SHEPHERD HOOK 0.8
- f. SHEPHERD HOOK 1.0
- g. SHEPHERD HOOK 1.0 MOD.
- h. RENAL DOUBLE CURVE




- a. STRAIGHT SELECTIVE
- b. AMIR MOTARJEME®
- c. AMIR MOTAJEME® CANE
- d. REUTER
- e. MIKAELSSON




- a. COBRA 1
- b. COBRA 2
- c. COBRA 3



- a. SHEPHERD FLUSH
- b. ULTRA BOLUS FLUSH
- c. MODIFIED BOLUS FLUSH
- d. STANDARD FLUSH
- e. STRAIGHT FLUSH
- f. MODIFIED HOOK FLUSH
- g. RBI
- h. RIM
- i. MINI RIM
- j. MINI RIM ST
- k. UTERINE ARTERY



- a. MODIFIED HOOK SELECTIVE 1
- b. MODIFIED HOOK SELECTIVE 2
- c. MODIFIED HOOK SELECTIVE 3



IMPRESS® DIAGNOSTIC PERIPHERAL CATHETERS

The Impress line of diagnostic peripheral catheters has been designed to deliver radiopaque media to selected sites in the vascular system. The tips are tapered to minimise vessel trauma, the accordion hubs are designed to minimise kinking, and the proprietary braid pattern shafts provide maximum flexibility without compromising steerability.

IMPRESS® HEADHUNTER CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length
------------------	-------------	-------------------	------------------------	--------


HEADHUNTER 1

410035HH1	4	.040"/1.02 mm	.035" (0.89 mm)	100 cm
410038HH1	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm
412535HH1	4	.040"/1.02 mm	.035" (0.89 mm)	125 cm
510035HH1	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm
510038HH1	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm
512535HH1	5	.046"/1.17 mm	.035" (0.89 mm)	125 cm
512538HH1	5	.046"/1.17 mm	.038" (0.97 mm)	125 cm



HEADHUNTER 3

410035HH3	4	.040"/1.02 mm	.035" (0.89 mm)	100 cm
410038HH3	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm
510035HH3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm
510038HH3	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm



IMPRESS® NEWTON CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length
------------------	-------------	-------------------	------------------------	--------

NEWTON 1

410038N1	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm
510038N1	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm

NEWTON 2


410038N2	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm
510038N2	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm

NEWTON 3

410038N3	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm
510038N3	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm

NEWTON 4

410038N4	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm
510038N4	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm




IMPRESS® BENTSON CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length
BENTSON 1				
46535B1	4	.040"/1.02 mm	.035" (0.89 mm)	65 cm
410035B1	4	.040"/1.02 mm	.035" (0.89 mm)	100 cm
410038B1	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm
56535B1	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm
510035B1	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm
510038B1	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm



BENTSON 2

410035B2	4	.040"/1.02 mm	.035" (0.89 mm)	100 cm
510035B2	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm



IMPRESS® MANI CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length
410035MANI	4	.040"/1.02 mm	.035" (0.89 mm)	100 cm
410038MANI	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm
510035MANI	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm
510038MANI	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm




IMPRESS® VERTEBRAL CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length
46538VER	4	.040"/1.02 mm	.038" (0.97 mm)	65 cm
410035VER	4	.040"/1.02 mm	.035" (0.89 mm)	100 cm
410038VER	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm
412535VER	4	.040"/1.02 mm	.035" (0.89 mm)	125 cm
56538VER	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm
510035VER	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm
510038VER	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm
512535VER	5	.046"/1.17 mm	.035" (0.89 mm)	125 cm



IMPRESS® MODIFIED CEREBRAL CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length
410035MC	4	.040"/1.02 mm	.035" (0.89 mm)	100 cm
410038MC	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm
510035MC	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm
510038MC	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm



IMPRESS® BERENSTEIN CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length
46535BER	4	.040"/1.02 mm	.035" (0.89 mm)	65 cm
46538BER	4	.040"/1.02 mm	.038" (0.97 mm)	65 cm
410035BER	4	.040"/1.02 mm	.035" (0.89 mm)	100 cm
410038BER	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm
412538BER	4	.040"/1.02 mm	.038" (0.97 mm)	125 cm
56535BER	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm
56538BER	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm
510035BER	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm
510038BER	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm
512538BER	5	.046"/1.17 mm	.038" (0.97 mm)	125 cm




IMPRESS® MW2 CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length
44038MW2	4	.040"/1.02 mm	.038" (0.97 mm)	40 cm
56538MW2	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm



IMPRESS® SIMMONS CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length
SIMMONS 1				
46535SIM1	4	.040"/1.02 mm	.035" (0.89 mm)	65 cm
46538SIM1	4	.040"/1.02 mm	.038" (0.97 mm)	65 cm
410035SIM1	4	.040"/1.02 mm	.035" (0.89 mm)	100 cm
410038SIM1	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm
56535SIM1	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm
56538SIM1	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm
510035SIM1	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm
510038SIM1	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm



SIMMONS 2

46538SIM2	4	.040"/1.02 mm	.038" (0.97 mm)	65 cm
410035SIM2	4	.040"/1.02 mm	.035" (0.89 mm)	100 cm
410038SIM2	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm
56538SIM2	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm
510035SIM2	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm
510038SIM2	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm



MODIFIED SIMMONS

48035MODS	4	.040"/1.02 mm	.035" (0.89 mm)	80 cm
48038MODS	4	.040"/1.02 mm	.038" (0.97 mm)	80 cm



IMPRESS® COBRA CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports
------------------	-------------	-------------------	------------------------	--------	-----------------

COBRA 1 (SIDE HOLES)

465352CB1	4	.040"/1.02 mm	.035" (0.89 mm)	65 cm	2
465382CB1	4	.040"/1.02 mm	.038" (0.97 mm)	65 cm	2
565352CB1	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	2
565382CB1	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	2

COBRA 1 (WITHOUT SIDE HOLES)

46538CB1	4	.040"/1.02 mm	.038" (0.97 mm)	65 cm	-
410038CB1	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm	-
56538CB1	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	-
510038CB1	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	-

COBRA 2 (SIDE HOLES)

465352CB2	4	.040"/1.02 mm	.035" (0.89 mm)	65 cm	2
465382CB2	4	.040"/1.02 mm	.038" (0.97 mm)	65 cm	2
565352CB2	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	2
565382CB2	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	2

COBRA 2 (WITHOUT SIDE HOLES)

46535CB2	4	.040"/1.02 mm	.035" (0.89 mm)	65 cm	-
46538CB2	4	.040"/1.02 mm	.038" (0.97 mm)	65 cm	-
410038CB2	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm	-
56538CB2	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	-
510038CB2	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	-

COBRA 1



COBRA 2



IMPRESS® HOOK CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports
------------------	-------------	-------------------	------------------------	--------	-----------------

HOOK 0.8

465352HK0-8	4	.040"/1.02 mm	.035" (0.89 mm)	65 cm	2
565352HK0-8	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	2

HOOK 0.8



HOOK 1.0

465352HK1-0	4	.040"/1.02 mm	.035" (0.89 mm)	65 cm	2
565352HK1-0	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	2

HOOK 1.0



IMPRESS® SHEPHERD HOOK CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports
------------------	-------------	-------------------	------------------------	--------	-----------------

SHEPHERD HOOK 0.8

465352SH0-8	4	.040"/1.02 mm	.035" (0.89 mm)	65 cm	2
565352SH0-8	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	2

SHEPHERD HOOK 0.8



SHEPHERD HOOK 1.0

465352SH1-0	4	.040"/1.02 mm	.035" (0.89 mm)	65 cm	2
565352SH1-0	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	2

SHEPHERD HOOK 1.0



IMPRESS® RENAL DOUBLE CURVE CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports
------------------	-------------	-------------------	------------------------	--------	-----------------

46535RC	4	.040"/1.02 mm	.035" (0.89 mm)	65 cm	-
465352RC	4	.040"/1.02 mm	.035" (0.89 mm)	65 cm	2
565352RC	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	2

RENAL DOUBLE CURVE



IMPRESS® AMIR MOTARJEME® CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports
------------------	-------------	-------------------	------------------------	--------	-----------------

AMIR MOTARJEME

46535MOT	4	.040"/1.02 mm	.035" (0.89 mm)	65 cm	-
410035MOT	4	.040"/1.02 mm	.035" (0.89 mm)	100 cm	-
56535MOT	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	-
510035MOT	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	-

MOTARJEME



AMIR MOTARJEME CANE

46535MOTC	4	.040"/1.02 mm	.035" (0.89 mm)	65 cm	-
410035MOTC	4	.040"/1.02 mm	.035" (0.89 mm)	100 cm	-

MOTARJEME CANE



IMPRESS® REUTER CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length
48038RE	4	.040"/1.02 mm	.038" (0.97 mm)	80 cm
58035RE	5	.046"/1.17 mm	.035" (0.89 mm)	80 cm
58038RE	5	.046"/1.17 mm	.038" (0.97 mm)	80 cm



REUTER

IMPRESS® MIKAELSSON CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length
48038MIK	4	.040"/1.02 mm	.038" (0.97 mm)	80 cm
58038MIK	5	.046"/1.17 mm	.038" (0.97 mm)	80 cm



MIKAELSSON

IMPRESS® KA2 CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length
44038KA2	4	.040"/1.02 mm	.038" (0.97 mm)	40 cm
46538KA2	4	.040"/1.02 mm	.038" (0.97 mm)	65 cm
54038KA2	5	.046"/1.17 mm	.038" (0.97 mm)	40 cm
56538KA2	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm



KA2

IMPRESS® UTERINE ARTERY CATHETER

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length
59038UAC	5	.046"/1.17 mm	.038" (0.97 mm)	90 cm



UTERINE ARTERY

IMPRESS® HOCKEY STICK CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Tip Length
410038HS	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm	4.5 cm
510038HS	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	4.5 cm



HOCKEY STICK

IMPRESS® MODIFIED HOOK SELECTIVE CATHETERS


Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length
MODIFIED HOOK 1				
48035MHK1	4	.040"/1.02 mm	.035" (0.89 mm)	80 cm
58035MHK1	5	.046"/1.17 mm	.035" (0.89 mm)	80 cm



MODIFIED HOOK 2				
48035MHK2	4	.040"/1.02 mm	.035" (0.89 mm)	80 cm
58035MHK2	5	.046"/1.17 mm	.035" (0.89 mm)	80 cm



MODIFIED HOOK 3				
48035MHK3	4	.040"/1.02 mm	.035" (0.89 mm)	80 cm
58035MHK3	5	.046"/1.17 mm	.035" (0.89 mm)	80 cm



IMPRESS® RBI (RENAL, BIFURCATION, ILIAC) CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports
465354RBI	4	.040"/1.02 mm	.035" (0.89 mm)	65 cm	4
565354RBI	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	4
56535RBI	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	-
510035RBI	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	-
5100354RBI	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	4



IMPRESS® RIM - MINI RIM SHORT-TIP CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports
46535RIM	4	.040"/1.02 mm	.035" (0.89 mm)	65 cm	-
56535RIM	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	-
56535MRIMST	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	-




IMPRESS® STRAIGHT SELECTIVE CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length
44035STS	4	.040"/1.02 mm	.035" (0.89 mm)	40 cm
46538STS	4	.040"/1.02 mm	.038" (0.97 mm)	65 cm
411035STS	4	.040"/1.02 mm	.035" (0.89 mm)	110 cm
411038STS	4	.040"/1.02 mm	.038" (0.97 mm)	110 cm
54038STS	5	.046"/1.17 mm	.038" (0.97 mm)	40 cm
56538STS	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm
511035STS	5	.046"/1.17 mm	.035" (0.89 mm)	110 cm
511038STS	5	.046"/1.17 mm	.038" (0.97 mm)	110 cm



IMPRESS® OSBORN CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length
410035OSB1	4	.040"/1.02 mm	0.035" (.089 mm)	100 cm
510035OSB1	5	.046"/1.17 mm	0.035" (.089 mm)	100 cm




IMPRESS® MULTIPURPOSE A1 CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length
410038MPA1	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm
412538MPA1	4	.040"/1.02 mm	.038" (0.97 mm)	125 cm
56538MPA1	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm
510038MPA1	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm
512538MPA1	5	.046"/1.17 mm	.038" (0.97 mm)	125 cm




IMPRESS® 30 CM CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Description
43035MODCB1	4	.040"/1.02 mm	.035" (0.89 mm)	30 cm	-	
43035BER	4	.040"/1.02 mm	.035" (0.89 mm)	30 cm	-	
43035KA2	4	.040"/1.02 mm	.035" (0.89 mm)	30 cm	-	
43035STR	4	.040"/1.02 mm	.035" (0.89 mm)	30 cm	-	
43035RIM	4	.040"/1.02 mm	.035" (0.89 mm)	30 cm	-	
4303510PIG-NB	4	.040"/1.02 mm	.035" (0.89 mm)	30 cm	10	non-braided
53035MODCB1	5	.046"/1.17 mm	.035" (0.89 mm)	30 cm	-	
53035BER	5	.046"/1.17 mm	.035" (0.89 mm)	30 cm	-	
53035KA2	5	.046"/1.17 mm	.035" (0.89 mm)	30 cm	-	
53035STR	5	.046"/1.17 mm	.035" (0.89 mm)	30 cm	-	
53035RIM	5	.046"/1.17 mm	.035" (0.89 mm)	30 cm	-	
5303510PIG-NB	5	.049"/1.25 mm	.035" (0.89 mm)	30 cm	10	non-braided




IMPRESS® SHEPHERD FLUSH CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Description
465356SHP-NB	4	.040"/1.02 mm	.035" (0.89 mm)	65 cm	6	non-braided
490356SHP-NB	4	.040"/1.02 mm	.035" (0.89 mm)	90 cm	6	non-braided
565356SHP-NB	5	.049"/1.25 mm	.035" (0.89 mm)	65 cm	6	non-braided



IMPRESS® ULTRA BOLUS FLUSH CATHETERS


Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Description
4653510UBF-NB	4	.040"/1.02 mm	.035" (0.89 mm)	65 cm	10	non-braided
4903510UBF-NB	4	.040"/1.02 mm	.035" (0.89 mm)	90 cm	10	non-braided
5653510UBF-NB	5	.049"/1.25 mm	.035" (0.89 mm)	65 cm	10	non-braided
5903510UBF-NB	5	.049"/1.25 mm	.035" (0.89 mm)	90 cm	10	non-braided



ULTRA BOLUS FLUSH

IMPRESS® MODIFIED BOLUS FLUSH CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Description
4653810MUBF-NB	4	.040"/1.02 mm	.038" (0.97 mm)	65 cm	10	non-braided
5653810MUBF-NB	5	.049"/1.25 mm	.038" (0.97 mm)	65 cm	10	non-braided



MODIFIED BOLUS FLUSH

IMPRESS® PIGTAIL FLUSH CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Description
4303510PIG-NB	4	.040"/1.02 mm	.035" (0.89 mm)	30 cm	10	non-braided
4653510PIG-NB	4	.040"/1.02 mm	.035" (0.89 mm)	65 cm	10	non-braided
4903510PIG-NB	4	.040"/1.02 mm	.035" (0.89 mm)	90 cm	10	non-braided
41103510PIG-NB	4	.040"/1.02 mm	.035" (0.89 mm)	110 cm	10	non-braided
5303510PIG-NB	5	.049"/1.25 mm	.035" (0.89 mm)	30 cm	10	non-braided
5653510PIG-NB	5	.049"/1.25 mm	.035" (0.89 mm)	65 cm	10	non-braided
5903510PIG-NB	5	.049"/1.25 mm	.035" (0.89 mm)	90 cm	10	non-braided
51103510PIG-NB	5	.049"/1.25 mm	.035" (0.89 mm)	110 cm	10	non-braided
5653810PIG-NB	5	.049"/1.25 mm	.038" (0.97 mm)	65 cm	10	non-braided
5903810PIG-NB	5	.049"/1.25 mm	.038" (0.97 mm)	90 cm	10	non-braided
51103810PIG-NB	5	.049"/1.25 mm	.038" (0.97 mm)	110 cm	10	non-braided



PIGTAIL

IMPRESS® STRAIGHT FLUSH CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Description
4653510STF-NB	4	.040"/1.02 mm	.035" (0.89 mm)	65 cm	10	non-braided
4903510STF-NB	4	.040"/1.02 mm	.035" (0.89 mm)	90 cm	10	non-braided
5653510STF-NB	5	.049"/1.25 mm	.035" (0.89 mm)	65 cm	10	non-braided
5903510STF-NB	5	.049"/1.25 mm	.035" (0.89 mm)	90 cm	10	non-braided



STRAIGHT

IMPRESS® MODIFIED HOOK FLUSH CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Description
465356MHK-NB	4	.040"/1.02 mm	.035" (0.89 mm)	65 cm	6	non-braided
490356MHK-NB	4	.040"/1.02 mm	.035" (0.89 mm)	90 cm	6	non-braided



MODIFIED HOOK FLUSH

IMPRESS® HYDROPHILIC DIAGNOSTIC PERIPHERAL CATHETERS

Impress hydrophilic catheters feature a hydrophilic coating designed to reduce friction for smoother navigation through tortuous anatomy. The hydrophilic coating is resistant to abrasion and is designed to maintain lubricity throughout the procedure. The tips are tapered to minimise vessel trauma, the accordion hubs are designed to minimise kinking, and the stainless steel braided shaft construction provides maximum flexibility without compromising steerability.

IMPRESS® HYDROPHILIC BENTSON


Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Coating Length
------------------	-------------	-------------------	------------------------	--------	----------------

BENTSON 1

410038B1-H	4	.40"/1.02 mm	.038" (0.97 mm)	100 cm	40 cm
510038B1-H	5	.40"/1.02 mm	.038" (0.97 mm)	100 cm	40 cm

BENTSON 2


410038B2-H	4	.40"/1.02 mm	.038" (0.97 mm)	100 cm	40 cm
510038B2-H	5	.40"/1.02 mm	.038" (0.97 mm)	100 cm	40 cm



IMPRESS® HYDROPHILIC BERENSTEIN

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Coating Length
------------------	-------------	-------------------	------------------------	--------	----------------

46538BER-H	4	.046"/1.17 mm	.038" (0.97 mm)	65 cm	25 cm
410038BER-H	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm	40 cm
412538BER-H	4	.040"/1.02 mm	.038" (0.97 mm)	125 cm	40 cm
56538BER-H	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	25 cm
510038BER-H	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	40 cm
512538BER-H	5	.038"/0.97 mm	.038" (0.97 mm)	125 cm	40 cm



IMPRESS® HYDROPHILIC COBRA


Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Coating Length
------------------	-------------	-------------------	------------------------	--------	----------------

COBRA 1

46538CB1-H	4	.040"/1.02 mm	.038" (0.97 mm)	65 cm	25 cm
56538CB1-H	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	25 cm

COBRA 2

46538CB2-H	4	.040"/1.02 mm	.038" (0.97 mm)	65 cm	25 cm
410038CB2-H	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm	40 cm
56538CB2-H	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	25 cm
510038CB2-H	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	40 cm



IMPRESS® HYDROPHILIC HEADHUNTER 1

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Coating Length
410038HH1-H	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm	40 cm
412538HH1-H	4	.040"/1.02 mm	.038" (0.97 mm)	125 cm	40 cm
510038HH1-H	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	40 cm
512538HH1-H	5	.046"/1.17 mm	.038" (0.97 mm)	125 cm	40 cm



IMPRESS® HYDROPHILIC KA2

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Coating Length
44038KA2-H	4	.040"/1.02 mm	.038" (0.97 mm)	40 cm	25 cm
46538KA2-H	4	.040"/1.02 mm	.038" (0.97 mm)	65 cm	25 cm
54038KA2-H	5	.046"/1.17 mm	.038" (0.97 mm)	40 cm	25 cm
56538KA2-H	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	25 cm



IMPRESS® HYDROPHILIC MANI

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Coating Length
410038MANI-H	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm	40 cm
510038MANI-H	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	40 cm



IMPRESS® HYDROPHILIC MULTIPURPOSE

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Coating Length
410038MPA1-H	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm	40 cm
510038MPA1-H	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	40 cm




IMPRESS® HYDROPHILIC RBI

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Coating Length	
56538RBI-H	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	25 cm	 RBI


IMPRESS® HYDROPHILIC SIMMONS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Coating Length	
------------------	-------------	-------------------	------------------------	--------	----------------	--

SIMMONS 1

410038SIM1-H	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm	40 cm	 SIMMONS 1
510038SIM1-H	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	40 cm	


SIMMONS 2

410038SIM2-H	4	.040" (1.02 mm)	.038" (0.97 mm)	100 cm	40 cm	 SIMMONS 2
510038SIM2-H	5	.046" (1.17 mm)	.038" (0.97 mm)	100 cm	40 cm	


IMPRESS® HYDROPHILIC MODIFIED HOOK

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Coating Length	
------------------	-------------	-------------------	------------------------	--------	----------------	--


MODIFIED HOOK 1

48038MHK1-H	4	.040"/1.02 mm	.038" (0.97 mm)	80 cm	40 cm	 MODIFIED HOOK 1
58038MHK1-H	5	.046"/1.17 mm	.038" (0.97 mm)	80 cm	40 cm	


MODIFIED HOOK 2

48038MHK2-H	4	.040"/1.02 mm	.038" (0.97 mm)	80 cm	40 cm	 MODIFIED HOOK 2
58038MHK2-H	5	.046"/1.17 mm	.038" (0.97 mm)	80 cm	40 cm	


MODIFIED HOOK 3

48038MHK3-H	4	.040"/1.02 mm	.038" (0.97 mm)	80 cm	40 cm	 MODIFIED HOOK 3
58038MHK3-H	5	.046"/1.17 mm	.038" (0.97 mm)	80 cm	40 cm	

IMPRESS® HYDROPHILIC STRAIGHT SELECTIVE

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Coating Length	
46538STS-H	4	.040"/1.02 mm	.038" (0.97 mm)	65 cm	25 cm	 <p>STRAIGHT SELECTIVE</p>
410038STS-H	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm	40 cm	
412538STS-H	4	.040"/1.02 mm	.038" (0.97 mm)	125 cm	40 cm	
56538STS-H	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	25 cm	
510038STS-H	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	40 cm	

IMPRESS® HYDROPHILIC VERTEBRAL





Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Coating Length	
46538VER-H	4	.040"/1.02 mm	.038" (0.97 mm)	65 cm	25 cm	 <p>VERTEBRAL</p>
410038VER-H	4	.040"/1.02 mm	.038" (0.97 mm)	100 cm	40 cm	
56538VER-H	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	25 cm	
510038VER-H	5	.046"/1.17 mm.	.038" (0.97 mm)	100 cm	40 cm	

PERFORMA® DIAGNOSTIC PERIPHERAL CATHETERS

Merit Medical's Performa diagnostic peripheral catheters are designed to provide superior torque, increased shaft strength for better pushability, and large inner lumens for improved flow rates.

PERFORMA® HEADHUNTER CATHETERS


Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Description
HEADHUNTER 1					
7715-A0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	
7715-10	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	bumper tip
5575-A3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	
5575-13	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	
HEADHUNTER 3					
7715-B0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	
7715-30	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	bumper tip
5575-B3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	
5575-23	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	
HEADHUNTER 5					
7715-C0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	
7715-50	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	bumper tip
5575-C3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	
5575-33	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	
HEADHUNTER 6					
7715-D0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	
7715-60	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	bumper tip
5575-D3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	
5575-43	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	



ANGIOGRAPHY

PERFORMA® NEWTON CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Description
NEWTON 1					
7722-A0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	
7722-10	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	bumper tip
5577-A3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	
5577-13	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	
NEWTON 2					
7722-B0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	
7722-20	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	bumper tip
5577-B3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	
5577-23	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	
NEWTON 3					
7722-C0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	
7722-30	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	bumper tip
5577-C3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	
5577-33	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	
NEWTON 4					
7722-D0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	
7722-40	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	bumper tip
5577-D3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	
5577-43	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	
NEWTON 5					
7722-E0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	
7722-50	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	bumper tip
5577-E3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	
5577-53	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	




PERFORMA® BENTSON CATHETERS


Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Description
BENTSON 1					
7721-A0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	
7721-10	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	bumper tip
5578-A3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	
5578-13	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	



BENTSON 2					
7721-B0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	
7721-20	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	bumper tip
5578-B3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	
5578-23	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	



BENTSON 3					
7721-C0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	
7721-30	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	bumper tip
5578-C3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	
5578-33	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	



PERFORMA® MANI CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Description
7716-A0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	
7716-10	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	bumper tip
5579-A3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	
5579-13	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	




PERFORMA® VERTEBRAL CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Description
7717-A0	4	.042"/1.07 mm	.035" (0.89 mm)	65 cm	
7717-10	4	.042"/1.07 mm	.038" (0.97 mm)	65 cm	bumper tip
7717-C0	4	.042"/1.07 mm	.035" (0.89 mm)	80 cm	
7717-30	4	.042"/1.07 mm	.038" (0.97 mm)	80 cm	bumper tip
7717-B0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	
7717-20	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	bumper tip
5583-B3	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	
5583-23	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	
5583-A3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	
5583-13	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	



PERFORMA® MODIFIED CEREBRAL CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Description
7723-A0	4	.042"/1.07 mm	.035" (0.89 mm)	90 cm	
7723-10	4	.042"/1.07 mm	.038" (0.97 mm)	90 cm	bumper tip
5588-A3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	
5588-13	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	



PERFORMA® BERENSTEIN CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Description
7718-C0	4	.042"/1.07 mm	.035" (0.89 mm)	40 cm	
7718-30	4	.042"/1.07 mm	.038" (0.97 mm)	40 cm	bumper tip
7718-B0	4	.042"/1.07 mm	.035" (0.89 mm)	65 cm	
7718-20	4	.042"/1.07 mm	.038" (0.97 mm)	65 cm	bumper tip
7718-A0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	
7718-10	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	bumper tip
5592-C3	5	.046"/1.17 mm	.035" (0.89 mm)	40 cm	
5592-33	5	.046"/1.17 mm	.038" (0.97 mm)	40 cm	
5592-B3	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	
5592-23	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	
5592-A3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	
5592-13	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	



PERFORMA® STRAIGHT SELECTIVE CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Description
7748-A0	4	.042"/1.07 mm	.035" (0.89 mm)	110 cm	
7748-13	4	.042"/1.07 mm	.038" (0.97 mm)	110 cm	bumper tip



PERFORMA® MW2 CATHETERS

MW2

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Description
7745-10	4	.042"/1.07 mm	.038" (0.97 mm)	40 cm	2	bumper tip
5590-13	5	.046"/1.17 mm	.038" (0.97 mm)	40 cm	-	
5590-23	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	-	
5590-33	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	-	

MODIFIED MW2

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Description
7745-A0	4	.042"/1.07 mm	.035" (0.89 mm)	40 cm	2	
7745-B0	4	.042"/1.07 mm	.035" (0.89 mm)	65 cm	2	
7745-C0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	2	



MW2



MODIFIED MW2

PERFORMA® OSBORN CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Description
------------------	-------------	-------------------	------------------------	--------	-------------

OSBORN 1

7746-A0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	
7746-10	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	bumper tip

OSBORN 2

7746-B0	4	.042"/1.07 mm	.035" (0.89 mm)	65 cm	
7746-20	4	.042"/1.07 mm	.038" (0.97 mm)	65 cm	bumper tip
7746-C0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	
7746-30	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	bumper tip



OSBORN 1



OSBORN 2

PERFORMA® HOOK CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Description
------------------	-------------	-------------------	------------------------	--------	-----------------	-------------

HOOK 0.8

7711-A0	4	.042"/1.07 mm	.035" (0.89 mm)	65 cm	2	
7711-10	4	.042"/1.07 mm	.038" (0.97 mm)	65 cm	2	bumper tip
5585-A3	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	2	
5585-13	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	2	

HOOK 1.0

7711-B0	4	.042"/1.07 mm	.035" (0.89 mm)	65 cm	2	
7711-20	4	.042"/1.07 mm	.038" (0.97 mm)	65 cm	2	bumper tip
5585-B3	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	2	
5585-23	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	2	




HOOK 0.8



HOOK 1.0


PERFORMA® MODIFIED HOOK SELECTIVE CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length
MODIFIED HOOK 1				
5585-C3	5	.046"/1.17 mm	.035" (0.89 mm)	80 cm
MODIFIED HOOK 2				
5585-D3	5	.046"/1.17 mm	.035" (0.89 mm)	80 cm
MODIFIED HOOK 3				
5585-E3	5	.046"/1.17 mm	.035" (0.89 mm)	80 cm



PERFORMA® SIMMONS CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Description
SIMMONS 1					
7720-A0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	
7720-10	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	bumper tip
5576-F3	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	
5576-A3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	
5576-13	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	
SIMMONS 2					
7720-B0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	
7720-20	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	bumper tip
5576-G3	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	
5576-B3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	
5576-23	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	
SIMMONS 3					
7720-C0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	
7720-30	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	bumper tip
5576-C3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	
5576-33	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	
SIMMONS 4					
7720-D0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	
7720-40	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	bumper tip
5576-D3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	
5576-43	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	
MODIFIED SIMMONS					
7720-53	4	.042"/1.07 mm	.038" (0.97 mm)	80 cm	bumper tip
MODIFIED SIMMONS 1.5					
7720-F0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	
7720-E0	4	.042"/1.07 mm	.035" (0.89 mm)	80 cm	
7720-70	4	.042"/1.07 mm	.038" (0.97 mm)	110 cm	bumper tip
7720-80	4	.042"/1.07 mm	.038" (0.97 mm)	120 cm	bumper tip
5576-E3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	




PERFORMA® COBRA CATHETERS

ANGIOGRAPHY

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Description
------------------	-------------	-------------------	------------------------	--------	-----------------	-------------


COBRA 1

7710-A0	4	.042"/1.07 mm	.035" (0.89 mm)	65 cm	2	
7710-10	4	.042"/1.07 mm	.038" (0.97 mm)	65 cm	2	bumper tip
7710-D0	4	.042"/1.07 mm	.035" (0.89 mm)	65 cm	-	
7710-40	4	.042"/1.07 mm	.038" (0.97 mm)	65 cm	-	bumper tip
5584-A3	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	2	
5584-D3	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	-	
5584-13	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	-	
5584-43	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	2	




COBRA 2

7710-B0	4	.042"/1.07 mm	.035" (0.89 mm)	65 cm	2	
7710-E0	4	.042"/1.07 mm	.035" (0.89 mm)	65 cm	-	
7710-20	4	.042"/1.07 mm	.038" (0.97 mm)	65 cm	2	bumper tip
7710-50	4	.042"/1.07 mm	.038" (0.97 mm)	65 cm	-	bumper tip
5584-B3	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	2	
5584-E3	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	-	
5584-23	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	-	
5584-53	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	2	



COBRA 3

7710-C0	4	.042"/1.07 mm	.035" (0.89 mm)	65 cm	2	
7710-F0	4	.042"/1.07 mm	.035" (0.89 mm)	65 cm	-	
7710-30	4	.042"/1.07 mm	.038" (0.97 mm)	65 cm	2	bumper tip
7710-60	4	.042"/1.07 mm	.038" (0.97 mm)	65 cm	-	bumper tip
5584-C3	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	2	
5584-F3	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	-	
5584-33	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	-	
5584-63	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	2	



PERFORMA® SHEPHERD HOOK CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Description
------------------	-------------	-------------------	------------------------	--------	-----------------	-------------

SHEPHERD HOOK 0.8

7712-A0	4	.042"/1.07 mm	.035" (0.89 mm)	65 cm	2	
7712-10	4	.042"/1.07 mm	.038" (0.97 mm)	65 cm	2	bumper tip
5586-A3	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	2	
5586-13	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	2	

SHEPHERD HOOK 0.8



SHEPHERD HOOK 1.0

7712-B0	4	.042"/1.07 mm	.035" (0.89 mm)	65 cm	2	
7712-20	4	.042"/1.07 mm	.038" (0.97 mm)	65 cm	2	bumper tip
5586-B3	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	2	
5586-23	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	2	

SHEPHERD HOOK 1.0



SHEPHERD HOOK 1.0 MODIFIED

7712-C0	4	.042"/1.07 mm	.035" (0.89 mm)	65 cm	2	
7712-D0	4	.042"/1.07 mm	.035" (0.89 mm)	65 cm	-	
7712-30	4	.042"/1.07 mm	.038" (0.97 mm)	65 cm	2	bumper tip
7712-40	4	.042"/1.07 mm	.038" (0.97 mm)	65 cm	-	bumper tip
5586-43	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	2	

SHEPHERD HOOK 1.0 MOD.



PERFORMA® RENAL DOUBLE CURVE CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Description
------------------	-------------	-------------------	------------------------	--------	-----------------	-------------

7713-A0	4	.042"/1.07 mm	.035" (0.89 mm)	65 cm	2	
7713-B0	4	.042"/1.07 mm	.035" (0.89 mm)	65 cm	-	
7713-10	4	.042"/1.07 mm	.038" (0.97 mm)	65 cm	2	bumper tip
7713-20	4	.042"/1.07 mm	.038" (0.97 mm)	65 cm	-	bumper tip
5587-A3	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	2	
5587-13	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	-	
5587-23	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	2	

RENAL DOUBLE CURVE



PERFORMA® HOCKEY STICK CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Tip Length	# of Side Ports	Description
------------------	-------------	-------------------	------------------------	--------	------------	-----------------	-------------

7752-10	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	4.5 cm	-	bumper tip
7543-13	5	.046"/1.17 mm	.038" (0.97 mm)	100 cm	4.5 cm	-	
7541-13	6	.054"/1.37 mm	.038" (0.97 mm)	100 cm	4.5 cm	-	

HOCKEY STICK



PERFORMA® AMIR MOTARJEME® CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Description
------------------	-------------	-------------------	------------------------	--------	-----------------	-------------

AMIR MOTARJEME


7747-C0	4	.042"/1.07 mm	.035" (0.89 mm)	65 cm	-	
7747-D0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	-	
7747-E0	4	.042"/1.07 mm	.035" (0.89 mm)	65 cm	2	
7747-F0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	2	
5591-A3	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	-	
5591-B3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	-	
5591-C3	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	2	
5591-D3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	2	



MOTARJEME

AMIR MOTARJEME CANE

7747-B0	4	.042"/1.07 mm	.035" (0.89 mm)	65 cm	-	
7747-20	4	.042"/1.07 mm	.038" (0.97 mm)	65 cm	-	bumper tip
7747-A0	4	.042"/1.07 mm	.035" (0.89 mm)	100 cm	-	
7747-10	4	.042"/1.07 mm	.038" (0.97 mm)	100 cm	-	bumper tip
5591-E3	5	.046"/1.17 mm	.035" (0.89 mm)	65 cm	-	
5591-F3	5	.046"/1.17 mm	.035" (0.89 mm)	100 cm	-	




MOTARJEME CANE

PERFORMA® REUTER CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Description
------------------	-------------	-------------------	------------------------	--------	-----------------	-------------

7714-A0	4	.042"/1.07 mm	.035" (0.89 mm)	80 cm	2	
7714-10	4	.042"/1.07 mm	.038" (0.97 mm)	80 cm	2	bumper tip
7714-B0	4	.042"/1.07 mm	.035" (0.89 mm)	80 cm	-	
7714-20	4	.042"/1.07 mm	.038" (0.97 mm)	80 cm	-	bumper tip
5589-A3	5	.046"/1.17 mm	.035" (0.89 mm)	80 cm	-	
5589-B3	5	.046"/1.17 mm	.035" (0.89 mm)	80 cm	2	
5589-13	5	.046"/1.17 mm	.038" (0.97 mm)	80 cm	-	
5589-23	5	.046"/1.17 mm	.038" (0.97 mm)	80 cm	2	



REUTER

PERFORMA® MIKAELSSON CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Description
7719-A0	4	.042"/1.07 mm	.035" (0.89 mm)	80 cm	
7719-10	4	.042"/1.07 mm	.038" (0.97 mm)	80 cm	bumper tip
5594-13	5	.046"/1.17 mm	.038" (0.97 mm)	80 cm	



MIKAELSSON

PERFORMA® KA CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	Description
5582-23	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm	
5582-33	5	.052"/1.32 mm	.038" (0.97 mm)	65 cm	high-flow, non-braided



KA

PERFORMA® KA2 CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length
5582-13	5	.046"/1.17 mm	.038" (0.97 mm)	40 cm
5582-43	5	.046"/1.17 mm	.038" (0.97 mm)	65 cm
7514-13	6	.054"/1.37 mm	.038" (0.97 mm)	40 cm
7514-23	6	.054"/1.37 mm	.038" (0.97 mm)	65 cm



KA2

PERFORMA® DVS A1 CATHETERS (10 CM TAPER TIP)

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length
5595-13	5	.046"/1.12 mm	.035" (0.89 mm)	65 cm
5595-33	5	.046"/1.12 mm	.035" (0.89 mm)	100 cm
5595-23	5	.046"/1.12 mm	.038" (0.97 mm)	65 cm
5595-43	5	.046"/1.12 mm	.038" (0.97 mm)	100 cm



DVS A1

PERFORMA® DVS A2 CATHETERS (15 CM TAPER TIP)


Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length
5595-A3	5	.046"/1.12 mm	.035" (0.89 mm)	65 cm
5595-C3	5	.046"/1.12 mm	.035" (0.89 mm)	100 cm
5595-B3	5	.046"/1.12 mm	.038" (0.97 mm)	65 cm
5595-D3	5	.046"/1.12 mm	.038" (0.97 mm)	100 cm



DVS A2

PERFORMA® UHF SHEPHERD FLUSH CATHETERS


Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Description
7602-A1	5	.052"/1.32 mm	.035" (0.89 mm)	65 cm	6	non-braided
7602-B1	5	.052"/1.32 mm	.035" (0.89 mm)	90 cm	6	non-braided
7602-C1	5	.052"/1.32 mm	.038" (0.97 mm)	65 cm	6	non-braided
7602-D1	5	.052"/1.32 mm	.038" (0.97 mm)	90 cm	6	non-braided



SHEPHERD FLUSH

PERFORMA® ULTRA BOLUS FLUSH CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Description
7674-11	5	.052"/1.32 mm	.035" (0.89 mm)	65 cm	10	non-braided
7674-21	5	.052"/1.32 mm	.035" (0.89 mm)	90 cm	10	non-braided
7674-31	5	.052"/1.32 mm	.038" (0.97 mm)	65 cm	10	non-braided
7674-41	5	.052"/1.32 mm	.038" (0.97 mm)	90 cm	10	non-braided



ULTRA BOLUS FLUSH

PERFORMA® ULTRA HIGH FLOW PIGTAIL CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Description
7602-41	5	.052"/1.32 mm	.035" (0.89 mm)	65 cm	10	non-braided
7602-51	5	.052"/1.32 mm	.035" (0.89 mm)	90 cm	10	non-braided
7602-61	5	.052"/1.32 mm	.035" (0.89 mm)	110 cm	10	non-braided
7602-71	5	.052"/1.32 mm	.038" (0.97 mm)	65 cm	10	non-braided
7602-81	5	.052"/1.32 mm	.038" (0.97 mm)	90 cm	10	non-braided
7602-91	5	.052"/1.32 mm	.038" (0.97 mm)	110 cm	10	non-braided



PIGTAIL

PERFORMA® PIGTAIL FLUSH CATHETER

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Description
7749-A0	4	.042"/1.07 mm	.035" (0.89 mm)	80 cm	10	braided



PIGTAIL

PERFORMA® STRAIGHT FLUSH CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Description
7673-11	5	.052"/1.32 mm	.035" (0.89 mm)	65 cm	10	non-braided
7673-21	5	.052"/1.32 mm	.035" (0.89 mm)	90 cm	10	non-braided
7673-31	5	.052"/1.32 mm	.035" (0.89 mm)	65 cm	-	non-braided
7673-41	5	.052"/1.32 mm	.038" (0.97 mm)	65 cm	10	non-braided
7673-51	5	.052"/1.32 mm	.038" (0.97 mm)	90 cm	10	non-braided
7673-61	5	.052"/1.32 mm	.038" (0.97 mm)	65 cm	-	non-braided

STRAIGHT



PERFORMA® MODIFIED HOOK FLUSH CATHETERS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Description
7603-A1	5	.052"/1.32 mm	.035" (0.89 mm)	65 cm	6	non-braided
7603-B1	5	.052"/1.32 mm	.035" (0.89 mm)	90 cm	6	non-braided
7603-C1	5	.052"/1.32 mm	.038" (0.97 mm)	65 cm	6	non-braided
7603-D1	5	.052"/1.32 mm	.038" (0.97 mm)	90 cm	6	non-braided

MODIFIED HOOK FLUSH





PERFORMA® VESSEL-SIZING CATHETERS

ALL SIDE PORTS BETWEEN TIP AND FIRST BAND

Vessel-Sizing Catheters (marker catheters) are designed to provide accurate measurement of the patient’s anatomy and will take the guesswork out of selecting appropriate interventional devices. Making the correct sizing decisions save time and produces better patient outcomes.

ANGIOGRAPHY

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Bands and Spacing	Description
2 BAND							
7742-2M65SH	4	.041"/1.04 mm	.035" (0.89 mm)	65 cm	5	2 (2 cm apart)	Pigtail
7742-2M100SH	4	.041"/1.04 mm	.035" (0.89 mm)	100 cm	5	2 (2 cm apart)	Pigtail
7744-2M65SH	4	.041"/1.04 mm	.035" (0.89 mm)	65 cm	5	2 (2 cm apart)	Straight
7744-2M100SH	4	.041"/1.04 mm	.035" (0.89 mm)	100 cm	5	2 (2 cm apart)	Straight
7743-2M65SH	4	.041"/1.04 mm	.035" (0.89 mm)	65 cm	5	2 (2 cm apart)	Ultra Bolus
7743-2M100SH	4	.041"/1.04 mm	.035" (0.89 mm)	100 cm	5	2 (2 cm apart)	Ultra Bolus
7755-2M65	4	.041"/1.04 mm	.035" (0.89 mm)	65 cm	6	2 (2 cm apart)	Mod. Hook
7755-2M100	4	.041"/1.04 mm	.035" (0.89 mm)	100 cm	6	2 (2 cm apart)	Mod. Hook
7602-2M65	5	.052"/1.32 mm	.035" (0.89 mm)	65 cm	10	2 (2 cm apart)	Pigtail
7602-2M	5	.052"/1.32 mm	.035" (0.89 mm)	100 cm	10	2 (2 cm apart)	Pigtail
7674-2M65	5	.052"/1.32 mm	.035" (0.89 mm)	65 cm	10	2 (2 cm apart)	Ultra Bolus
7674-2M	5	.052"/1.32 mm	.035" (0.89 mm)	100 cm	10	2 (2 cm apart)	Ultra Bolus
7673-2M65	5	.052"/1.32 mm	.035" (0.89 mm)	65 cm	10	2 (2 cm apart)	Straight
7673-2M	5	.052"/1.32 mm	.035" (0.89 mm)	100 cm	10	2 (2 cm apart)	Straight
7603-2M65	5	.052"/1.32 mm	.035" (0.89 mm)	65 cm	6	2 (2 cm apart)	Mod. Hook
7603-2M100	5	.052"/1.32 mm	.035" (0.89 mm)	100 cm	6	2 (2 cm apart)	Mod. Hook



11 BAND

7602-11M65	5	.052"/1.32 mm	.035" (0.89 mm)	65 cm	10	11 (1 cm apart)	Pigtail
7602-11M	5	.052"/1.32 mm	.035" (0.89 mm)	100 cm	10	11 (1 cm apart)	Pigtail
7673-11M65SH	5	.052"/1.32 mm	.035" (0.89 mm)	65 cm	10	11 (1 cm apart)	Straight
7673-11M100SH	5	.052"/1.32 mm	.035" (0.89 mm)	100 cm	10	11 (1 cm apart)	Straight
7603-11M65	5	.052"/1.32 mm	.035" (0.89 mm)	65 cm	6	11 (1 cm apart)	Mod. Hook
7603-11M100	5	.052"/1.32 mm	.035" (0.89 mm)	100 cm	6	11 (1 cm apart)	Mod. Hook



20 BAND

7742-20M65SH	4	.041"/1.04 mm	.035" (0.89 mm)	65 cm	5	20 (1 cm apart)	Pigtail
7742-20M100SH	4	.041"/1.04 mm	.035" (0.89 mm)	100 cm	5	20 (1 cm apart)	Pigtail
7602-20M65	5	.052"/1.32 mm	.035" (0.89 mm)	65 cm	10	20 (1 cm apart)	Pigtail
7602-20M	5	.052"/1.32 mm	.035" (0.89 mm)	100 cm	10	20 (1 cm apart)	Pigtail
7602-20M65SH	5	.052"/1.32 mm	.035" (0.89 mm)	65 cm	5	20 (1 cm apart)	Pigtail
7602-20M100SH	5	.052"/1.32 mm	.035" (0.89 mm)	100 cm	5	20 (1 cm apart)	Pigtail
7673-20M65SH	5	.052"/1.32 mm	.035" (0.89 mm)	65 cm	10	20 (1 cm apart)	Straight
7673-20M100SH	5	.052"/1.32 mm	.035" (0.89 mm)	100 cm	10	20 (1 cm apart)	Straight
7603-20M65SH	5	.052"/1.32 mm	.035" (0.89 mm)	65 cm	6	20 (1 cm apart)	Mod. Hook
7603-20M100SH	5	.052"/1.32 mm	.035" (0.89 mm)	100 cm	6	20 (1 cm apart)	Mod. Hook




SH suffix identifies marker bands that are 16 mm closer to the distal end.

PERFORMA® VESSEL-SIZING CATHETERS

ALL SIDE PORTS BETWEEN BANDS

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Bands and Spacing	Description
2 BAND**							
56535AVS2PIG	5	.052"/1.32 mm	.035" (0.89 mm)	65 cm	10	2 (2 cm apart)	Pigtail
510035AVS2PIG	5	.052"/1.32 mm	.035" (0.89 mm)	100 cm	10	2 (2 cm apart)	Pigtail
56535AVS2STR	5	.052"/1.32 mm	.035" (0.89 mm)	65 cm	10	2 (2 cm apart)	Straight
510035AVS2STR	5	.052"/1.32 mm	.035" (0.89 mm)	100 cm	10	2 (2 cm apart)	Straight
56535AVS2MHK	5	.052"/1.32 mm	.035" (0.89 mm)	65 cm	6	2 (2 cm apart)	Mod. Hook
510035AVS2MHK	5	.052"/1.32 mm	.035" (0.89 mm)	100 cm	6	2 (2 cm apart)	Mod. Hook



2 BAND STRAIGHT



2 BAND ULTRA BOLUS




2 BAND MOD. HOOK FLUSH



Standard Configuration

Platinum marker bands are placed at the distal portion of the catheter to determine accurate vessel sizing.



SH Configuration

Platinum marker bands are positioned closer to the pigtail curve allowing bands to be visualised closer to the vessel.



AVS Configuration

Platinum marker bands have side ports placed between bands to allow quick contrast filling of the vessel while visualising the marker bands.

Catalogue Number	French Size	Internal Diameter	Recommended Guide Wire	Length	# of Side Ports	Bands and Spacing	Description
20 BAND**							
46535AVS20MHK	4	.041"/1.04 mm	.035" (0.89 mm)	65 cm	6	20 (1 cm apart)	Mod. Hook
56535AVS20PIG	5	.052"/1.32 mm	.035" (0.89 mm)	65 cm	10	20 (1 cm apart)	Pigtail
510035AVS20PIG	5	.052"/1.32 mm	.035" (0.89 mm)	100 cm	10	20 (1 cm apart)	Pigtail
56535AVS20MHK	5	.052"/1.32 mm	.035" (0.89 mm)	65 cm	6	20 (1 cm apart)	Mod. Hook
510035AVS20MHK	5	.052"/1.32 mm	.035" (0.89 mm)	100 cm	6	20 (1 cm apart)	Mod. Hook



20 BAND PIGTAIL

** Platinum marker bands have side ports placed between bands to allow quick contrast filling of the vessel while visualising the marker bands.

HIGH PRESSURE FLOW CONTROL SWITCH

The Merit Medical High Pressure Flow Control Switch gives you the confidence you need when stopping and starting flow during critical procedures. This feature-rich, integrated, one-handed, single-channel device is designed with clinician and patient safety in mind. Packaged 24 units per box.



Catalogue Number U1SHC

CT TRANSFER SETS


Merit Medical's line of CT Transfer Sets are interchangeable for use with either contrast or saline solutions. The CT Transfer Set features a wide base for insertion support when placing into the contrast bottle. Available with a large bore piercing spike and a large wing, swabbable valve. The saline transfer set has a longer, sharper spike tip for easier insertion into IV solution bags and a blue safety stripe for identification.



Catalogue Number

TS120

Contrast spike with 20" (50 cm) small bore tubing and HALKEY-ROBERTS® brand swabbable valve. Packaged 50/box, 4 boxes/case.



Catalogue Number

TS220


Saline spike with 20" (50 cm) small bore blue stripe tubing and HALKEY-ROBERTS® brand swabbable valve. Packaged 20/box, 4 boxes/case.



Catalogue Number

TS520L

Saline spike with 20" (50 cm) large bore tubing and HALKEY-ROBERTS® brand swabbable valve. Packaged 20/box, 4 boxes/case.



Catalogue Number

K08-01961B

Vented, valved spike with 20" (50 cm) small bore tubing and nexus Luer. Packaged 25/box, 4 boxes/case.



Catalogue Number

K08-02195A

Vented, valved spike with 20" (50 cm) small bore blue stripe tubing and nexus Luer. Packaged 25/box, 4 boxes/case.



HALKEY-ROBERTS® is a registered trademark of Halkey-Roberts Corporation.

MERIT MAESTRO® MICROCATHERETERS

The Merit Maestro is a multipurpose microcatheter designed for use in peripheral and coronary vasculature. It features a hydrophilic coating for smooth traction through tortuous anatomy, a radiopaque marker for easy identification of the distal tip, and an exclusive swan neck design to help seat the catheter in the vessel.



Catalogue Number	Proximal	Distal	Length cm	Tip	Minimum Inner Diameter	Recommended Guide Wire	Recommended Guiding Catheter
2.8F TAPERED TO 2.4F							
28MC24110ST	2.8	2.4	110	Straight	.020" (0.52 mm)	.018" (0.46 mm)	.040" (1.02 mm)
28MC24130ST	2.8	2.4	130	Straight	.020" (0.52 mm)	.018" (0.46 mm)	.040" (1.02 mm)
28MC24150ST	2.8	2.4	150	Straight	.020" (0.52 mm)	.018" (0.46 mm)	.040" (1.02 mm)
28MC2411045	2.8	2.4	110	45 degree	.020" (0.52 mm)	.018" (0.46 mm)	.040" (1.02 mm)
28MC2413045	2.8	2.4	130	45 degree	.020" (0.52 mm)	.018" (0.46 mm)	.040" (1.02 mm)
28MC2415045	2.8	2.4	150	45 degree	.020" (0.52 mm)	.018" (0.46 mm)	.040" (1.02 mm)
28MC24110SN	2.8	2.4	110	Swan Neck	.020" (0.52 mm)	.018" (0.46 mm)	.040" (1.02 mm)
28MC24130SN	2.8	2.4	130	Swan Neck	.020" (0.52 mm)	.018" (0.46 mm)	.040" (1.02 mm)
28MC24150SN	2.8	2.4	150	Swan Neck	.020" (0.52 mm)	.018" (0.46 mm)	.040" (1.02 mm)

2.8F NO TAPER

28MC28110ST	2.8	2.8	110	Straight	.024" (0.62 mm)	.021" (0.53 mm)	.040" (1.02 mm)
28MC28130ST	2.8	2.8	130	Straight	.024" (0.62 mm)	.021" (0.53 mm)	.040" (1.02 mm)
28MC28150ST	2.8	2.8	150	Straight	.024" (0.62 mm)	.021" (0.53 mm)	.040" (1.02 mm)
28MC2811045	2.8	2.8	110	45 degree	.024" (0.62 mm)	.021" (0.53 mm)	.040" (1.02 mm)
28MC2813045	2.8	2.8	130	45 degree	.024" (0.62 mm)	.021" (0.53 mm)	.040" (1.02 mm)
28MC2815045	2.8	2.8	150	45 degree	.024" (0.62 mm)	.021" (0.53 mm)	.040" (1.02 mm)
28MC28110SN	2.8	2.8	110	Swan Neck	.024" (0.62 mm)	.021" (0.53 mm)	.040" (1.02 mm)
28MC28130SN	2.8	2.8	130	Swan Neck	.024" (0.62 mm)	.021" (0.53 mm)	.040" (1.02 mm)
28MC28150SN	2.8	2.8	150	Swan Neck	.024" (0.62 mm)	.021" (0.53 mm)	.040" (1.02 mm)

2.9F NO TAPER

29MC29110ST	2.9	2.9	110	Straight	.027" (0.68 mm)	.021" (0.53 mm)	.042" (1.07 mm)
29MC29130ST	2.9	2.9	130	Straight	.027" (0.68 mm)	.021" (0.53 mm)	.042" (1.07 mm)
29MC29150ST	2.9	2.9	150	Straight	.027" (0.68 mm)	.021" (0.53 mm)	.042" (1.07 mm)
29MC2911045	2.9	2.9	110	45 degree	.027" (0.68 mm)	.021" (0.53 mm)	.042" (1.07 mm)
29MC2913045	2.9	2.9	130	45 degree	.027" (0.68 mm)	.021" (0.53 mm)	.042" (1.07 mm)
29MC2915045	2.9	2.9	150	45 degree	.027" (0.68 mm)	.021" (0.53 mm)	.042" (1.07 mm)
29MC29110SN	2.9	2.9	110	Swan Neck	.027" (0.68 mm)	.021" (0.53 mm)	.042" (1.07 mm)
29MC29130SN	2.9	2.9	130	Swan Neck	.027" (0.68 mm)	.021" (0.53 mm)	.042" (1.07 mm)
29MC29150SN	2.9	2.9	150	Swan Neck	.027" (0.68 mm)	.021" (0.53 mm)	.042" (1.07 mm)


All Merit Maestro Microcatheters are rated to 800 psi infusion dynamic pressure. Microcatheters sold individually; 1 per box.

MEDALLION® SYRINGES

Merit Medical's Medallion syringes come in seven colour options with custom printing, or a frosted field for writing information directly on the syringe barrel. The Medallion Syringe features a clear polycarbonate barrel, coloured plunger, and fixed male Luer. Packaged 25 units per box, 4 boxes per case.




MSS011-LG



Plunger Colours

White	MSS011	MSS031	MSS061	MSS111	MSS121	MSS131	MSS160
Light Blue	MSS011-LB	MSS031-LB	MSS061-LB	MSS111-LB	MSS121-LB	MSS131-LB	
Yellow	MSS011-Y	MSS031-Y	MSS061-Y	MSS111-Y	MSS121-Y	MSS131-Y	
Red	MSS011-R	MSS031-R	MSS061-R	MSS111-R	MSS121-R	MSS131-R	
Light Green	MSS011-LG	MSS031-LG	MSS061-LG	MSS111-LG	MSS121-LG	MSS131-LG	
Dark Green	MSS011-DG	MSS031-DG	MSS061-DG	MSS111-DG	MSS121-DG	MSS131-DG	
Purple	MSS011-PR	MSS031-PR	MSS061-PR	MSS111-PR	MSS121-PR	MSS131-PR	

MEDALLION SYRINGES with clear polycarbonate barrel, coloured plunger, black medication printing, and fixed male Luer. Packaged 25 units per box, 4 boxes per case.



Plunger Colours


	NITRO (10 mL)	LIDOCAINE (10 mL)	HEPARIN (10 mL)	CONTRAST (10 mL)	SALINE (10 mL)
Light Blue	K01-60085	K01-60081	K01-60082	K01-60086	K01-60083
Yellow	K01-60101	K01-60097	K01-60098	K01-60102	K01-60099
Red	K01-60109	K01-60105	K01-60106	K01-60110	
White	K01-60077			K01-60078	K01-60075
Light Green	K01-60093	K01-60089	K01-60090	K01-60094	
Dark Green				K01-00786	

Plunger Colours

	FENTANYL (3 mL)	XYLOCAINE (10 mL)	VERSED (3 mL)	FLUSH (10 mL)
Light Blue		K01-00917	K01-05185	K01-05093
Red		K01-00953		
Light Green	K01-03243			

MEDALLION® SWORD HANDLE SYRINGES

Medallion Sword Handle Syringes feature sword handles, clear polycarbonate barrels, coloured plungers, and fixed male Luers. Packaged 25 units per box, 4 boxes per case.



Plunger Colours	10 mL	20 mL
White	MSSW11	MSSW21
Light Blue	MSSW11-LB	MSSW21-LB
Yellow	MSSW11-Y	MSSW21-Y
Red	MSSW11-R	MSSW21-R
Light Green	MSSW11-LG	MSSW21-LG
Dark Green	MSSW11-DG	MSSW21-DG
Purple	MSSW11-PR	MSSW21-PR



MSSW11-Y

MEDALLION® FROSTED FIELD SYRINGES & PEN

The Medallion Frosted Field Syringe allows doctors or clinicians to write additional pertinent information on the barrel. The syringe features a clear polycarbonate barrel, white plunger, white frosted field, and fixed male Luer. Packaged 25 units per box, 4 boxes per case.



Catalogue Number	1 mL	3 mL	6 mL	10 mL	20 mL	30 mL	60 mL
	MSSF011	MSSF031	MSSF061	MSSF111	MSSF121	MSSF131	MSSF160

PEN for labeling Medallion Frosted Field Syringes.



Catalogue Number K13-00476



Packaged 25 units per box, 4 boxes per case.

CUSTOM MEDALLION® SYRINGES

Merit Medical’s signature Medallion Syringes can be custom ordered to your specifications. Medallion Syringes come with a clear polycarbonate barrel, coloured plunger, customised medication printing, and Luer options including fixed male or slip Luer styles.

Select one option from each category and identify your custom medication to be printed on the barrel. Your Merit Representative will provide you with the corresponding catalogue number.

Custom Medallion Syringes are packaged 25 units per box, 4 boxes per case.



SIZE OPTIONS 1 mL 3 mL 6 mL 10 mL 20 mL 30 mL 60 mL

PLUNGER COLOUR OPTIONS:

- White
- Light Blue
- Light Green
- Dark Green
- Yellow
- Red
- Purple

*60 mL plunger is only available in white, light blue or red.

PRINT COLOUR OPTIONS:

- Black
- White

LUER OPTIONS:

- Fixed Male Luer
- Slip Luer



VACLOK® VACUUM PRESSURE SYRINGES

The VacLok vacuum pressure syringe comes with a clear polycarbonate barrel, white plunger, four- or six-position locking mechanism and fixed male Luer. Packaged 20 or 25 units per box, 4 boxes per case.



10 mL
four position
25/box

20 mL
four position
20/box

30 mL
six position
20/box

60 mL
six position
25/box

Catalogue Number

K01-07736

K01-07738

K01-07739

K01-07740