



*Merit Medical's Innovative Family
of Vascular Access Products*

A Great Introduction

Prelude®

SHEATH INTRODUCERS

The Prelude family of sheath introducers is a comprehensive product line offering excellent transitions, smooth insertion and more.

HEMOSTASIS VALVE - Minimizes blood leakage. Facilitates insertion and withdrawal of devices.

Catalog Number	French Size	Length (cm)	Guide Wire	Guide Wire Diameter	Wire Length	Hub Color
----------------	-------------	-------------	------------	---------------------	-------------	-----------

PRELUDE SHEATHS

PSI-4F-11-035	4F	11cm	Yes	0.035" (0.89mm)	50cm	red
PSI-5F-11-035	5F	11cm	Yes	0.035" (0.89mm)	50cm	gray
PSI-6F-11-035	6F	11cm	Yes	0.035" (0.89mm)	50cm	green
PSI-7F-11-035	7F	11cm	Yes	0.035" (0.89mm)	50cm	orange
PSI-8F-11-035	8F	11cm	Yes	0.035" (0.89mm)	50cm	blue

PSI-5F-11-038	5F	11cm	Yes	0.038" (0.97mm)	50cm	gray
PSI-6F-11-038	6F	11cm	Yes	0.038" (0.97mm)	50cm	green
PSI-7F-11-038	7F	11cm	Yes	0.038" (0.97mm)	50cm	orange
PSI-8F-11-038	8F	11cm	Yes	0.038" (0.97mm)	50cm	blue

PSI-4F-11	4F	11cm	No	-	-	red
PSI-5F-11	5F	11cm	No	-	-	gray
PSI-6F-11	6F	11cm	No	-	-	green
PSI-7F-11	7F	11cm	No	-	-	orange
PSI-8F-11	8F	11cm	No	-	-	blue

PSI-4F-23-035	4F	23cm	Yes	0.035" (0.89mm)	80cm	red
PSI-5F-23-035	5F	23cm	Yes	0.035" (0.89mm)	80cm	gray
PSI-6F-23-035	6F	23cm	Yes	0.035" (0.89mm)	80cm	green
PSI-7F-23-035	7F	23cm	Yes	0.035" (0.89mm)	80cm	orange
PSI-8F-23-035	8F	23cm	Yes	0.035" (0.89mm)	80cm	blue

PSI-5F-23-038	5F	23cm	Yes	0.038" (0.97mm)	80cm	gray
PSI-6F-23-038	6F	23cm	Yes	0.038" (0.97mm)	80cm	green
PSI-7F-23-038	7F	23cm	Yes	0.038" (0.97mm)	80cm	orange
PSI-8F-23-038	8F	23cm	Yes	0.038" (0.97mm)	80cm	blue

PSI-4F-23	4F	23cm	No	-	-	red
PSI-5F-23	5F	23cm	No	-	-	gray
PSI-6F-23	6F	23cm	No	-	-	green
PSI-7F-23	7F	23cm	No	-	-	orange
PSI-8F-23	8F	23cm	No	-	-	blue

Packaged in holster; 5 per box



Prelude^{PRO}

SHEATH INTRODUCERS

Polypropylene (PRO) sheath tubing offers increased resilience and additional support during procedure.

COLOR-CODED HUBS - Match standard French sizing for easy identification.

SNAP FIT DILATOR - Eliminates slipping during insertion. Easy to grasp surface with an audible snap.



Catalog Number	French Size	Length (cm)	Guide Wire	Guide Wire Diameter	Wire Length	Hub Color
----------------	-------------	-------------	------------	---------------------	-------------	-----------

PRELUDE PRO SHEATHS

PRO-4F-11-035	4F	11cm	Yes	0.035" (0.89mm)	50cm	red
PRO-5F-11-035	5F	11cm	Yes	0.035" (0.89mm)	50cm	gray
PRO-6F-11-035	6F	11cm	Yes	0.035" (0.89mm)	50cm	green
PRO-7F-11-035	7F	11cm	Yes	0.035" (0.89mm)	50cm	orange
PRO-8F-11-035	8F	11cm	Yes	0.035" (0.89mm)	50cm	blue

PRO-5F-11-038	5F	11cm	Yes	0.038" (0.97mm)	50cm	gray
PRO-6F-11-038	6F	11cm	Yes	0.038" (0.97mm)	50cm	green
PRO-7F-11-038	7F	11cm	Yes	0.038" (0.97mm)	50cm	orange
PRO-8F-11-038	8F	11cm	Yes	0.038" (0.97mm)	50cm	blue

PRO-4F-11	4F	11cm	No	-	-	red
PRO-5F-11	5F	11cm	No	-	-	gray
PRO-6F-11	6F	11cm	No	-	-	green
PRO-7F-11	7F	11cm	No	-	-	orange
PRO-8F-11	8F	11cm	No	-	-	blue

Packaged in holster; 5 per box



Prelude®

MARKER TIP SHEATH INTRODUCERS

Highly radiopaque marker provides increased visibility of the sheath tip when using fluoroscopy.



Catalog Number	French Size	Length (cm)	Guide Wire	Guide Wire Diameter	Wire Length	Hub Color
----------------	-------------	-------------	------------	---------------------	-------------	-----------

PRELUDE MT SHEATHS

▲ PRO-4F-11-035MT	4F	11cm	Yes	0.035" (0.89mm)	50cm	red
PRO-5F-11-035MT	5F	11cm	Yes	0.035" (0.89mm)	50cm	gray
PRO-6F-11-035MT	6F	11cm	Yes	0.035" (0.89mm)	50cm	green
PRO-7F-11-035MT	7F	11cm	Yes	0.035" (0.89mm)	50cm	orange
PRO-8F-11-035MT	8F	11cm	Yes	0.035" (0.89mm)	50cm	blue
PRO-5F-11-038MT	5F	11cm	Yes	0.038" (0.97mm)	50cm	gray
PRO-6F-11-038MT	6F	11cm	Yes	0.038" (0.97mm)	50cm	green
PRO-7F-11-038MT	7F	11cm	Yes	0.038" (0.97mm)	50cm	orange
PRO-8F-11-038MT	8F	11cm	Yes	0.038" (0.97mm)	50cm	blue
PRO-4F-11MT	4F	11cm	No	-	-	red
PRO-5F-11MT	5F	11cm	No	-	-	gray
PRO-6F-11MT	6F	11cm	No	-	-	green
PRO-7F-11MT	7F	11cm	No	-	-	orange
PRO-8F-11MT	8F	11cm	No	-	-	blue
PSI-4F-23-035MT	4F	23cm	Yes	0.035" (0.89mm)	80cm	red
PSI-5F-23-035MT	5F	23cm	Yes	0.035" (0.89mm)	80cm	gray
PSI-6F-23-035MT	6F	23cm	Yes	0.035" (0.89mm)	80cm	green
PSI-7F-23-035MT	7F	23cm	Yes	0.035" (0.89mm)	80cm	orange
PSI-8F-23-035MT	8F	23cm	Yes	0.035" (0.89mm)	80cm	blue
PSI-5F-23-038MT	5F	23cm	Yes	0.038" (0.97mm)	80cm	gray
PSI-6F-23-038MT	6F	23cm	Yes	0.038" (0.97mm)	80cm	green
PSI-7F-23-038MT	7F	23cm	Yes	0.038" (0.97mm)	80cm	orange
PSI-8F-23-038MT	8F	23cm	Yes	0.038" (0.97mm)	80cm	blue
PSI-4F-23MT	4F	23cm	No	-	-	red
PSI-5F-23MT	5F	23cm	No	-	-	gray
PSI-6F-23MT	6F	23cm	No	-	-	green
PSI-7F-23MT	7F	23cm	No	-	-	orange
PSI-8F-23MT	8F	23cm	No	-	-	blue

▲ All catalog numbers beginning with PRO are Prelude PRO tubing.

Packaged in holster; 5 per box

Prelude®


RADIAL (MINI ACCESS) SHEATH INTRODUCERS

Facilitates sheath placement via small needle (21g or 20g) and small guide wire (0.018"/0.46mm or 0.025"/0.64mm).

Product is packaged with needle, guide wire, sheath and dilator.

Excellent for radial, brachial, femoral and pediatric use.

EXCELLENT TRANSITIONS - Facilitates smooth entry into patient and vessel. Resists kinking during insertion. Dilators are available with 0.018" (0.46mm) and 0.025" (0.64mm) tip tapers.



Catalog Number	French Size	Length (cm)	Guide Wire Diameter & Length	Needle Gauge & Length	Sheath Hub Color
----------------	-------------	-------------	------------------------------	-----------------------	------------------

PRELUDE MINI ACCESS SHEATHS

PSI-4F-7-018	4F	7cm	0.018" (0.46mm) x 40cm SS/SS	21g x 4cm	red
PSI-5F-7-018	5F	7cm	0.018" (0.46mm) x 40cm SS/SS	21g x 4cm	gray
PSI-6F-7-018	6F	7cm	0.018" (0.46mm) x 40cm SS/SS	21g x 4cm	green
PSI-4F-11-018	4F	11cm	0.018" (0.46mm) x 40cm SS/SS	21g x 4cm	red
PSI-5F-11-018	5F	11cm	0.018" (0.46mm) x 40cm SS/SS	21g x 4cm	gray
PSI-6F-11-018	6F	11cm	0.018" (0.46mm) x 40cm SS/SS	21g x 4cm	green
PSI-4F-7-018NT4	4F	7cm	0.018" (0.46mm) x 40cm NT/PT	21g x 4cm	red
PSI-5F-7-018NT4	5F	7cm	0.018" (0.46mm) x 40cm NT/PT	21g x 4cm	gray
PSI-6F-7-018NT4	6F	7cm	0.018" (0.46mm) x 40cm NT/PT	21g x 4cm	green
PSI-4F-11-018NT4	4F	11cm	0.018" (0.46mm) x 40cm NT/PT	21g x 4cm	red
PSI-5F-11-018NT4	5F	11cm	0.018" (0.46mm) x 40cm NT/PT	21g x 4cm	gray
PSI-6F-11-018NT4	6F	11cm	0.018" (0.46mm) x 40cm NT/PT	21g x 4cm	green
▲ PRO-4F-11-018NT7	4F	11cm	0.018" (0.46mm) x 40cm NT/PT	21g x 7cm	red
PRO-5F-11-018NT7	5F	11cm	0.018" (0.46mm) x 40cm NT/PT	21g x 7cm	gray
PRO-6F-11-018NT7	6F	11cm	0.018" (0.46mm) x 40cm NT/PT	21g x 7cm	green
PSI-4F-23-018NT4	4F	23cm	0.018" (0.46mm) x 80cm NT/PT	21g x 4cm	red
PSI-5F-23-018NT4	5F	23cm	0.018" (0.46mm) x 80cm NT/PT	21g x 4cm	gray
PSI-6F-23-018NT4	6F	23cm	0.018" (0.46mm) x 80cm NT/PT	21g x 4cm	green
PSI-4F-23-018NT7	4F	23cm	0.018" (0.46mm) x 80cm NT/PT	21g x 7cm	red
PSI-5F-23-018NT7	5F	23cm	0.018" (0.46mm) x 80cm NT/PT	21g x 7cm	gray
PSI-6F-23-018NT7	6F	23cm	0.018" (0.46mm) x 80cm NT/PT	21g x 7cm	green

SS = Stainless Steel NT = Nitinol PT = Platinum tip

▲ All catalog numbers beginning with PRO are Prelude PRO tubing.

PSI-4F-7-025	4F	7cm	0.025" (0.64mm) x 50cm	20g x 4cm	red
PSI-5F-7-025	5F	7cm	0.025" (0.64mm) x 50cm	20g x 4cm	gray
PSI-6F-7-025	6F	7cm	0.025" (0.64mm) x 50cm	20g x 4cm	green
PSI-4F-11-025	4F	11cm	0.025" (0.64mm) x 50cm	20g x 4cm	red
PSI-5F-11-025	5F	11cm	0.025" (0.64mm) x 50cm	20g x 4cm	gray
PSI-6F-11-025	6F	11cm	0.025" (0.64mm) x 50cm	20g x 4cm	green

Guide wire is stainless steel; double ended (3mmJ and floppy).

Packaged in holster; 5 per box

MERIT **MAK**TM

MAK AND S-MAK MINI ACCESS KITS

*Excellent transition between dilator and introducer.
Dilator fits snugly into the introducer for smooth entry
into the vessel. Stiffened dilator (patent pending design)
provides additional support during placement through
scar tissue and/or difficult anatomy.*

*Kit configurations provide the necessary
components for quick access to the vessel.*

Catalog Number	French Size	Introducer Length	Guide Wire Diameter/Length	Guide Wire Core/Tip	Needle Gauge/Length	Qty/Box
----------------	-------------	-------------------	----------------------------	---------------------	---------------------	---------

STANDARD 7CM NEEDLE

MAK401	4F	10cm	*0.018" x 40cm	SS/SS	21G x 7cm	10
MAK501	5F	10cm	0.018" x 40cm	SS/SS	21G x 7cm	10
MAK401N	4F	10cm	0.018" x 40cm	NT/PT	21G x 7cm	10
MAK401NPD	4F	10cm	0.018" x 40cm	NT/PD	21G x 7cm	10
MAK501N	5F	10cm	0.018" x 40cm	NT/PT	21G x 7cm	10
MAK501NPD	5F	10cm	0.018" x 40cm	NT/PD	21G x 7cm	10
MAK501N-15	5F	15cm	0.018" x 60cm	NT/PT	21G x 7cm	10
MAK405	4F	10cm	0.018" x 40cm	SS/PT	21G x 7cm	10
MAK405PD	4F	10cm	0.018" x 40cm	SS/PD	21G x 7cm	10
MAK505	5F	10cm	0.018" x 40cm	SS/PT	21G x 7cm	10
MAK505PD	5F	10cm	0.018" x 40cm	SS/PD	21G x 7cm	10

ECHO ENHANCED 7CM NEEDLE

MAK401-E	4F	10cm	0.018" x 40cm	SS/SS	21G x 7cm	10
MAK501-E	5F	10cm	0.018" x 40cm	SS/SS	21G x 7cm	10
MAK401N-E	4F	10cm	0.018" x 40cm	NT/PT	21G x 7cm	10
MAK401NPD-E	4F	10cm	0.018" x 40cm	NT/PD	21G x 7cm	10
MAK501N-E	5F	10cm	0.018" x 40cm	NT/PT	21G x 7cm	10
MAK501NPD-E	5F	10cm	0.018" x 40cm	NT/PD	21G x 7cm	10
MAK405-E	4F	10cm	0.018" x 40cm	SS/PT	21G x 7cm	10
MAK505-E	5F	10cm	0.018" x 40cm	SS/PT	21G x 7cm	10

*0.018" (0.46mm) SS = Stainless Steel PT = Platinum Tip NT = Nitinol
PD = Palladium Tip

MERITMEDICALTM

Distribuido por:



CardioParts[®]
Tecnología con sentido humano

ventas@cardioparts.com • www.cardioparts.com
Tels. (55) 9183-3933 (81) 8357-0779 (81) 8357-1903

