

RAD BOARD® 2


RADIAL ARM BOARD

The Next Generation of Radial

Provides built-in radiation scatter protection



THINK RADIAL.
DISCOVER MERIT.™



Rad Board® 2

Radial Arm Board

OUR RADIAL ADVANTAGE

The Rad Board 2 is the latest addition to Merit's comprehensive radial portfolio.

- Designed with a sleek, smaller profile for clinical utility
- Radiation Scatter Protection
 - 23-25% neck height¹
 - 43-44% waist height¹
- Easy to use with bi-plane equipment
- Reversible and reusable
- Product dimensions: 28"/71 cm x 23"/58 cm



RELATED PRODUCTS



Rad Trac™

Allows easy placement of the Rad Board 2 when patient is on the table

Fits all cath lab/radiology tables

Reusable

Product dimensions:
26.5"/67 cm x 11"/28 cm



Rad Board Xtra™

Allows 90° perpendicular extension of arm during access

Provides an additional 11 inches (28 cm) of working space when placed next to Rad Board 2

Reversible and reusable

Product dimensions:
26.5"/67 cm x 11"/28 cm



Rad Rest®

Provides wrist and elbow support

Stabilizes and positions arm

Easily adjusts to accommodate any upper extremity length

Soft foam for increased patient comfort

Rad Board® 2

Radial Arm Board



The Rad Board and accessories set-up provides convenience, comfort and maximizes working space.



The Rad Trac fits under the patient and the mattress allowing for the Rad Board 2 to conveniently slide under the patient with ease.

The Rad Board can be used with or without the RadTrac or Rad Board Xtra.

Easy to position
under patient



NEXT GENERATION OF RADIAL

Rad Board® 2

Radial Arm Board

ORDERING INFORMATION

Catalog Number	Description	Order Unit
RB050	Rad Board 2 arm board	Sold individually, non-sterile
RBX	Rad Board Xtra	Sold individually, non-sterile
RT913	Rad Trac	Sold individually, non-sterile
RR811	Rad Rest	Individually packaged, 10 units per box, non-sterile

¹Data on file.




THINK RADIAL. DISCOVER MERIT.™

From set up to hemostasis, Merit Medical offers a complete radial product portfolio to help you achieve the procedural **EFFICIENCY** and **SAFETY** that is essential for delivering the best possible **CARE**.



SET-UP & ACCESS



ANGIOGRAPHY



INTERVENTION



HEMOSTASIS



For more information on our radial product portfolio, scan the QR code or visit merit.com/tr. For more information regarding our extensive radial education opportunities, visit thinkradial.com

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Merit.com

Distribuido por:

ventas@cardioparts.com • www.cardioparts.com
Tels. (55) 9183-3933 (81) 8357-0779 (81) 8357-1903