



**Quality Plus Customization
Equals Convenience**


TRAM™
INTEGRATED TRANSDUCER

YOUR BOTTOM LINE IS CONVENIENCE.

Enjoy the convenience of custom designing your own kit with our feature-rich products to support your unique clinical requirements. Merit Medical offers an array of high-value fluid management products and angiography accessories to help you save time and money and to keep pace with your procedural volumes.

Merit's new TRAM manifold with an integrated transducer conveniently located on the first port gives a single operator the ability to de-bubble, zero and re-zero within the sterile field. Increase comfort and reduce hand fatigue with our low-torque, easy-to-turn wide-spaced handles.





ORDERING INFORMATION

TRAM™ manifold with an integrated transducer on the port closest to the rotating adapter. Comfortable wide spaced finger grips. Packaged 10 units per box, 4 boxes per case.

	1 Port	2 Port	3 Port	4 Port
200 PSI (14 BAR)				
ON Handles	201LTPHN-R	202LTPHWN-R	203LTPHWN-R	204LTPHWN-R
OFF Handles	201LTPHF-R	202LTPHWF-R	203LTPHWF-R	204LTPHWF-R
500 PSI (34 BAR)				
ON Handles	501LTPHN-R	502LTPHWN-R	503LTPHWN-R	504LTPHWN-R
OFF Handles	501LTPHF-R	502LTPHWF-R	503LTPHWF-R	504LTPHWF-R

TRAM-P™ compensating manifold with an integrated transducer on the port closest to the rotating adapter with a pressure monitoring line off the compensation port to allow pressure readings at any level. Comfortable wide spaced finger grips. Packaged 10 units per box, 4 boxes per case.

	1 Port	2 Port	3 Port	4 Port
200 PSI (14 BAR)				
ON Handles	201LTCHN-R	202LTCHWN-R	203LTCHWN-R	204LTCHWN-R
OFF Handles	201LTCHF-R	202LTCHWF-R	203LTCHWF-R	204LTCHWF-R
500 PSI (34 BAR)				
ON Handles	501LTCHN-R	502LTCHWN-R	503LTCHWN-R	504LTCHWN-R
OFF Handles	501LTCHF-R	502LTCHWF-R	503LTCHWF-R	504LTCHWF-R

*48" Compensation Line and Clip

MANIFOLDS with right port orientation and rotating luer. Packaged 25 units per box, 4 boxes per case.

	2 Port Half Body	3 Port Half Body	4 Port Half Body	4 Port Block Body	3 Port Block Body	2 Port Block Body	1 Port Block Body
200 PSI (14 BAR)							
ON Handles	202HN-R	203HN-R	204HN-R	204BN-R	203BN-R	202BN-R	
Wide Port	202LHWN-R	203LHWN-R	204LHWN-R				
OFF Handles	202HF-R	203HF-R	204HF-R	204BF-R	203BF-R	202BF-R	
Wide Port	202LHWF-R	203LHWF-R	204LHWF-R				
360/OFF Handles				204BF-RT	203BF-RT		
500 PSI (34 BAR)							
ON Handles	502HN-R	503HN-R	504HN-R	504BN-R	503BN-R	502BN-R	501BN-R
Wide Port	502LHWN-R	503LHWN-R	504LHWN-R				
OFF Handles	502HF-R	503HF-R	504HF-R	504BF-R	503BF-R	502BF-R	501BF-R
Wide Port	502LHWF-R	503LHWF-R	504LHWF-R				

DEVOS™ MANIFOLDS Devos Manifolds with check relief valve, right port orientation and rotating luer. TRAM* and TRAM-P* have the added benefit of an integrated transducer on the port closest to the rotating adapter with comfortable wide spaced finger grips. Packaged 25 units per box, 4 boxes per case. *Packaged 10 units per box, 4 boxes per case.

	4 Port Half Body	3 Port Half Body	2 Port Half Body	4 Port TRAM	4 Port TRAM-P *
200 PSI (14 BAR)					
OFF Handles	2V4LHWF-R	2V3HF-R	2V2HF-R	2V4LTPHWF-R	2V4LTCHWF-R
360/OFF Handles	2V4LHWF-RT	2V3HF-RT	2V2HF-RT	2V4LTPHWF-RT	2V4LTCHWF-RT

*48" Compensation Line and Clip

Visit merit.com for
additional manifold
configurations

