

TEMNO ELITE™

Biopsy System

72%
LARGER
SAMPLES



BUILT FOR BETTER BIOPSY

The TEMNO Elite™ Biopsy System offers the reliability you expect with new features designed to provide superior soft tissue samples and total procedure control. An added patient-safety option minimizes the risk of complications, particularly during lung biopsies.² By listening to your feedback, Merit has modernized biopsy procedures with the TEMNO Elite—a device built for better biopsy.

Sharp Biopsy Needle

Designed to minimize tissue insertion force.

Sample Assist™ Feature

Reduces cumbersome sample retrieval with tweezers or a scalpel!

Adjustable Throw

(10 mm or 20 mm)
Provides precise sampling for a variety of procedural situations.

72%
LARGER
SAMPLES¹

Versatility for Your Procedures

| | | |
|-----|---|--|
| 14G | ● | Breast |
| 16G | ● | Liver, kidney, soft tissue masses |
| 18G | ● | Liver, kidney, thyroid, soft tissue masses |
| 20G | ● | Lung, thyroid, soft tissue masses |

Available in 6–20 cm lengths





SUPERIOR SAMPLES

Total Core™ Biopsy Technology
Average 72% larger than competitor side-notch, semi-automatic devices!¹



PROCEDURE CONTROL

Sample Assist
“Hands-free” feature retrieves full tissue samples with a push of a button!¹



PATIENT SAFETY

Valved Coaxial Introducer
Intended to minimize the risk of air entering the pleural space and fluid leakage²

Visit merit.com/TemnoElite for product demonstration and ordering information.

TEMNOELITE™ ORDERING INFORMATION

| Needle Gauge | Standalone Biopsy Needle | Biopsy Needle Size | Bundled Biopsy Needle with Valved Coaxial Introducer | Bundled Biopsy Needle with Non-valved Coaxial Introducer | Introducer Needle Size | Compatible Valved Coaxial Introducer | Compatible Non-valved Coaxial Introducer |
|--------------|--------------------------|--------------------|--|--|------------------------|--------------------------------------|--|
| 14G ● | TE1406 | 14G x 6 cm | - | - | | - | - |
| | TE1410 | 14G x 10 cm | TEVC1410 | TEC1410 | 13G x 5.3 cm | VC1410 | CX1410 |
| | TE1415 | 14G x 15 cm | TEVC1415 | TEC1415 | 13G x 10.3 cm | VC1415 | CX1415 |
| 16G ● | TE1606 | 16G x 6 cm | - | - | | - | - |
| | TE1610 | 16G x 10 cm | TEVC1610 | TEC1610 | 15G x 5.3 cm | VC1610 | CX1610 |
| | TE1615 | 16G x 15 cm | TEVC1615 | TEC1615 | 15G x 10.3 cm | VC1615 | CX1615 |
| 18G ● | TE1806 | 18G x 6 cm | - | - | | - | - |
| | TE1810 | 18G x 10 cm | TEVC1810 | TEC1810 | 17G x 5.3 cm | VC1810 | CX1810 |
| | TE1815 | 18G x 15 cm | TEVC1815 | TEC1815 | 17G x 10.3 cm | VC1815 | CX1815 |
| | TE1820 | 18G x 20 cm | TEVC1820 | TEC1820 | 17G x 15.3 cm | VC1820 | CX1820 |
| 20G ● | TE2006 | 20G x 6 cm | - | - | | - | - |
| | TE2010 | 20G x 10 cm | TEVC2010 | TEC2010 | 19G x 5.4 cm | VC2010 | CX2010 |
| | TE2015 | 20G x 15 cm | TEVC2015 | TEC2015 | 19G x 10.4 cm | VC2015 | CX2015 |
| | TE2020 | 20G x 20 cm | TEVC2020 | TEC2020 | 19G x 15.4 cm | VC2020 | CX2020 |

For single use only. Five devices per case.

REFERENCES: 1. Internal testing with 18G. Data on File. Marketing Claims for the TEMNO Elite Device (ECN147270, ECN152582). 2. Internal Testing. Data on File. Validated with data from Valved One-Step™ Centesis Catheter. See product brochure (401273001/C, ECN98478).

Built for Better Biopsy

The TEMNO Elite™ Biopsy System offers the reliability you expect with new features designed to provide superior soft tissue samples and total procedure control. An added patient-safety option minimizes the risk of complications, particularly during lung biopsies.² By listening to your feedback, Merit has modernized biopsy procedures with the TEMNO Elite—a device built for better biopsy.

Versatility for Your Procedures

| | |
|-----|--|
| 14G | ● Breast |
| 16G | ● Liver, kidney, soft tissue masses |
| 18G | ● Liver, kidney, thyroid, soft tissue masses |
| 20G | ● Lung, thyroid, soft tissue masses |

Available in 6–20 cm lengths

Adjustable Throw

(10 mm or 20 mm)
Provides precise sampling for a variety of procedural situations.

Sharp Biopsy Needle

Designed to minimize tissue insertion force.

Sample Assist™ Feature

Reduces cumbersome sample retrieval with tweezers or a scalpel!



SUPERIOR SAMPLES

Total Core™ Biopsy Technology
Average 72% larger than competitor side-notch, semi-automatic devices!



PROCEDURE CONTROL

Sample Assist
“Hands-free” feature retrieves full tissue samples with a push of a button!



PATIENT SAFETY

Valved Coaxial Introducer
Intended to minimize the risk of air entering the pleural space and fluid leakage?

Visit merit.com/TemnoElite for product demonstration and ordering information.

CREATE BIOPSY SOLUTIONS



BIOPSY PROCEDURE KITS

- K12T-04321** | Biopsy Procedure Basic Kit
- K12T-04322** | Biopsy Procedure Complete Kit

REFERENCES: 1. Internal testing with 18G. Data on File. Marketing Claims for the TEMNO Elite Device (ECN147270, ECN152582). 2. Internal Testing. Data on File. Validated with data from Valved One-Step™ Centesis Catheter. See product brochure (401273001/C, ECN98478).

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
Main: +1.801.253.1600
Customer Service: +1.801.208.4300

Distribuido por:

ventas@cardioparts.com • www.cardioparts.com
Tels. (55) 9183-3933 (81) 8357-0779 (81) 8357-1903